

Aflibercept ELISA Kit

Overview

Catalog No. DB941058

Specifications Aflibercept

Applications ELISA

Accession P15692

Stability and Storage When the kit was stored at the recommended temperature for 6 months, the signal intensity decreased by less than 20%. For unopened kits, if you want to prolong the storage time, please store the Standard, Detection A, Detection B and Microplate at - 20 °C, the rest reagents should be store at 4°C.

Detection method Colorimetric

Sample type Plasma, Serum

Assay type Quantitative

Sensitivity 40.12 ng/mL

Range 78.13 - 5, 000 ng/mL

Recovery 80-120%

Alternative name 862111-32-8

Shipping 2-8 °C

Note For Research Use Only.