

## Anti-Burosumab ELISA Kit

### Overview

---

<b>Catalog No.</b>	AV625018	
<b>Applications</b>	ELISA	
<b>Stability and Storage</b>	The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 10% prior to the expiration date under appropriate storage condition.	
<b>Detection method</b>	Colorimetric	
<b>Sample type</b>	Plasma, Serum	
<b>Assay type</b>	Quantitative	
<b>Recovery</b>	80%-120%	
<b>Product Usage Information</b>	Application	Dilution
	ELISA	1:2000-1:20000
<b>Note</b>	For Research Use Only.	