

Anti-Golimumab ELISA Kit

Overview

Catalog No.	AF879068
Specifications	Golimumab
Applications	ELISA
Accession	P01375
Stability and Storage	When the kit was stored at the recommended temperature for 6 months, the signal intensity decreased by less than 20%. For unopened kits, if you want to prolong the storage time, please store the Standard, Detection A, Detection B and Microplate at - 20°C, the rest reagents should be store at 4°C.
Detection method	Colorimetric
Sample type	Plasma, Serum
Assay type	Quantitative
Sensitivity	85.1 ng/mL
Range	156.25 - 10, 000 ng/mL
Recovery	80-120%
Alternative name	CNTO 148, 476181-74-5
Shipping	2-8 °C
Note	For Research Use Only.