

## Anti-MPXV M1R Human IgG ELISA Kit

### Overview

---

**Catalog No.** AK785028

**Specifications** M1R

**Applications** ELISA

**Stability and Storage** When the kit was stored at the recommended temperature for 6 months, the signal intensity decreased by less than 20%. For unopened kits, if you want to prolong the storage time, please store the Standard, Detection A, Detection B and Microplate at - 20 °C, the rest reagents should be store at 4°C.

**Detection method** Colorimetric

**Sample type** Plasma, Serum

**Assay type** Quantitative

**Sensitivity** 0.26 ng/mL

**Range** 0.31 - 20 ng/mL

**Recovery** 80-120%

**Shipping** 2-8 °C

**Note** For Research Use Only.