

Anti-Palivizumab ELISA Kit

Overview

Catalog No.	AK421018
Specifications	Palivizumab
Applications	ELISA
Accession	P03420
Stability and Storage	The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 10% prior to the expiration date under appropriate storage condition.
Detection method	Colorimetric
Sample type	Plasma, Serum
Assay type	Quantitative
Range	156.25 - 10, 000 ng/mL
Recovery	80-120%
Alternative name	MEDI-493, 188039-54-5
Shipping	2-8 °C
Note	For Research Use Only.