

Anti-SARS-CoV-2 RBD (BA.5) Human IgM ELISA Kit

Overview

| | |
|------------------------------|---|
| Catalog No. | AK565038 |
| Specifications | RBD |
| Applications | ELISA |
| Stability and Storage | The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 10% prior to the expiration date under appropriate storage condition. |
| Detection method | Colorimetric |
| Sample type | Plasma, Serum |
| Assay type | Quantitative |
| Sensitivity | 22.1 ng/mL |
| Range | 93.75 - 6, 000 ng/mL |
| Recovery | 80-120% |
| Shipping | 2-8 °C |
| Note | For Research Use Only. |