

Anti-SARS-CoV-2 Spike Protein Mouse IgA ELISA Kit

Overview

Catalog No.	VK474018
Specifications	S protein
Applications	ELISA
Stability and Storage	The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 10% prior to the expiration date under appropriate storage condition.
Detection method	Colorimetric
Sample type	Plasma, Serum
Assay type	Quantitative
Sensitivity	47.47 pg/mL
Range	125 - 8, 000 pg/mL
Recovery	80-120%
Shipping	2-8 °C
Note	For Research Use Only.