

Anti-Semaglutide ELISA Kit

Overview

Catalog No.	AP899018
Specifications	Semaglutide
Applications	ELISA
Stability and Storage	When the kit was stored at the recommended temperature for 6 months, the signal intensity decreased by less than 20%. For unopened kits, if you want to prolong the storage time, please store the Standard, Detection A, Detection B and Microplate at - 20 °C, the rest reagents should be store at 4°C.
Detection method	Colorimetric
Sample type	Plasma, Serum
Assay type	Quantitative
Sensitivity	6.49 ng/mL
Range	15.63 - 10, 000 ng/mL
Recovery	80-120%
Alternative name	CAS, 910463-68-2, Semaglutide, Rybelsus, Ozempic, NN9535, OG217SC, NNC 0113-0217
Shipping	2-8 °C
Note	For Research Use Only.