

## Anti-Tirzepatide (LY3298176) Human IgG ELISA Kit

### Overview

---

<b>Catalog No.</b>	AP899028
<b>Specifications</b>	Tirzepatide
<b>Applications</b>	ELISA
<b>Stability and Storage</b>	The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 10% prior to the expiration date under appropriate storage condition.
<b>Detection method</b>	Colorimetric
<b>Sample type</b>	Plasma, Serum
<b>Assay type</b>	Quantitative
<b>Sensitivity</b>	0.28 ng/mL
<b>Range</b>	0.31 - 20 ng/mL
<b>Recovery</b>	80-120%
<b>Alternative name</b>	Tirzepatide, LY3298176, 2023788-19-2
<b>Shipping</b>	2-8 °C
<b>Note</b>	For Research Use Only.