

Basiliximab ELISA Kit

Overview

Catalog No. DF996018

Specifications Basiliximab

Applications ELISA

Accession P01589

Stability and Storage When the kit was stored at the recommended temperature for 6 months, the signal intensity decreased by less than 20%. For unopened kits, if you want to prolong the storage time, please store the Standard, Detection A, Detection B and Microplate at - 20 °C, the rest reagents should be store at 4°C.

Detection method Colorimetric

Sample type Plasma, Serum

Assay type Quantitative

Sensitivity 13.31 ng/mL

Range 23.44 - 1, 500 ng/mL

Recovery 80-120%

Alternative name CHI621, CHI-621 / SDZ-CHI-621, 179045-86-4

Shipping 2-8 °C

Note For Research Use Only.