

Bimagrumab ELISA Kit

Overview

Catalog No. DB949018

Specifications Bimagrumab

Applications ELISA

Accession P27037 & Q13705

Stability and Storage When the kit was stored at the recommended temperature for 6 months, the signal intensity decreased by less than 20%. For unopened kits, if you want to prolong the storage time, please store the Standard, Detection A, Detection B and Microplate at - 20 °C, the rest reagents should be store at 4°C.

Detection method Colorimetric

Sample type Plasma, Serum

Assay type Quantitative

Sensitivity 97.28 ng/mL

Range 125 - 4, 000 ng/mL

Recovery 80-120%

Alternative name BYM338, 1356922-05-8

Shipping 2-8 °C

Note For Research Use Only.