

Brolucizumab ELISA Kit

Overview

Catalog No.	DB941028
Specifications	Brolucizumab
Applications	ELISA
Accession	P15692
Stability and Storage	The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 10% prior to the expiration date under appropriate storage condition.
Detection method	Colorimetric
Sample type	Plasma, Serum
Assay type	Quantitative
Range	1.5625 ng/mL – 100 ng/mL
Recovery	80-120%
Alternative name	ESBA-1008, ESBA1008, RTH258, 1531589-13-5
Shipping	2-8 °C
Note	For Research Use Only.