

Casirivimab ELISA kit

Overview

Catalog No. DK565028

Specifications Casirivimab

Applications ELISA

Accession P0DTC2

Stability and Storage When the kit was stored at the recommended temperature for 6 months, the signal intensity decreased by less than 20%. For unopened kits, if you want to prolong the storage time, please store the Standard, Detection A, Detection B and Microplate at - 20°C, the rest reagents should be store at 4°C.

Detection method Colorimetric

Sample type Plasma, Serum

Assay type Quantitative

Range 15.625 ng/mL - 1, 000 ng/mL

Recovery 80-120%

Alternative name REGN-10933, REGN10933, 2415933-42-3

Shipping 2-8 °C

Note For Research Use Only.