

# Certolizumab ELISA Kit

## Overview

---

**Catalog No.** DF879028

**Specifications** Certolizumab

**Applications** ELISA

**Accession** P01375

**Stability and Storage** When the kit was stored at the recommended temperature for 6 months, the signal intensity decreased by less than 20%. For unopened kits, if you want to prolong the storage time, please store the Standard, Detection A, Detection B and Microplate at - 20 °C, the rest reagents should be store at 4°C.

**Detection method** Colorimetric

**Sample type** Plasma, Serum

**Assay type** Quantitative

**Sensitivity** 329.66 ng/mL

**Range** 468.75 - 30, 000 ng/mL

**Recovery** 80-120%

**Alternative name** CDP870, PHA-738144, Certolizumab pegol, 428863-50-7

**Shipping** 2-8 °C

**Note** For Research Use Only.