

# Dupilumab ELISA Kit

## Overview

---

**Catalog No.** DB732018

**Specifications** Dupilumab

**Applications** ELISA

**Accession** P24394

**Stability and Storage** When the kit was stored at the recommended temperature for 6 months, the signal intensity decreased by less than 20%. For unopened kits, if you want to prolong the storage time, please store the Standard, Detection A, Detection B and Microplate at - 20°C, the rest reagents should be store at 4°C.

**Detection method** Colorimetric

**Sample type** Plasma, Serum

**Assay type** Quantitative

**Range** 78.125 ng/mL - 5, 000 ng/mL

**Recovery** 90-120%

**Alternative name** CD124, REGN668, SAR231893, 1190264-60-8

**Shipping** 2-8 °C

**Note** For Research Use Only.