

SARS-CoV-2 Surrogate Virus Neutralization Test Kit (B.1.617.1/E484Q, L452R)

Overview

Catalog No.	AK726018
Specifications	RBD (Kappa)
Applications	ELISA
Stability and Storage	The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 10% prior to the expiration date under appropriate storage condition.
Detection method	Colorimetric
Sample type	Plasma, Serum
Assay type	Quantitative
Sensitivity	78.1 ng/mL
Range	0.156 - 10 µg/mL
Recovery	80-120%
Alternative name	Spike glycoprotein, S glycoprotein, Spike protein S1, S, Receptor-Binding Domain, RBD, India-202010, Kappa, B.1.617.1, G/452R.V3, 21A/S:154K, Neutralization, SARS-CoV-2
Shipping	2-8 °C
Note	For Research Use Only.