

EasyBlot®


Make your perfect IP-WB blot easy

Your Expertise

Our Antibodies


Accelerated Discovery

EasyBlot®
Simple and Clear




Clear, efficiently eliminates denatured IgG masking of your target protein's signal
Easy, just switch from regular secondary antibodies to EasyBlot® secondary antibodies

Rabbit EasyBlot® with EasyBlocker (GTx225856-01)




IP analysis of p53 from 293T cell extracts using control rabbit IgG (GTx35035) or rabbit polyclonal anti-p53 antibody (GTx100629). WB analysis of the immunoprecipitated sample was performed using the same primary antibody (GTx100629).

Mouse EasyBlot® with EasyBlocker (GTx225857-01)



IP analysis of NUF2 from HeLa cell extracts using control mouse IgG (GTx35009) or mouse monoclonal anti-NUF2 antibody (GTx632184). WB analysis of the immunoprecipitated sample was performed using the same primary antibody (GTx632184).

EasyBlocker (GTx425858)



EasyBlocker, significantly reduce background generated by free Protein A/G! Western blot analysis of Protein A and Protein G samples. Membranes were blocked with 5% milk (A) or EasyBlocker (B). Membranes were probed with a Rabbit IgG followed by Rabbit EasyBlot.



EasyBlot Sampler Kit (GTX300111)

Citation Support


Cat. No.	Product Name	Package
GTX221666-01	EasyBlot anti Rabbit IgG (HRP)	10 µl
GTX221667-01	EasyBlot anti Mouse IgG (HRP)	10 µl
GTX425858	EasyBlocker	2 g

On any EasyBlot item


Citation Support Comparable Abs

Eliminate	Cat. No.	Product Name	Package	
Masking bands by denatured IgG	GTX221666-01	EasyBlot® anti-Rabbit IgG (HRP)	50 µl	250 µl
	GTX221667-01	EasyBlot® anti-Mouse IgG (HRP)		
	GTX628474-01	EasyBlot® anti-Rat IgG (HRP)		
	GTX628547-01	EasyBlot® anti-Goat IgG (HRP)		
	GTX628906-01	EasyBlot® anti-Sheep IgG (HRP)		
Masking bands by denatured IgG & Background from protein A/G	GTX225856-01	EasyBlot® anti-Rabbit IgG Kit (w/ EasyBlocker)	12 g	62 g
	GTX225857-01	EasyBlot® anti-Mouse IgG Kit (w/ EasyBlocker)		
	GTX228909-01	EasyBlot® anti-Rat IgG Kit (w/ EasyBlocker)		
	GTX228910-01	EasyBlot® anti-Goat IgG Kit (w/ EasyBlocker)		
	GTX228911-01	EasyBlot® anti-Sheep IgG Kit (w/ EasyBlocker)		
Background from protein A/G	GTX425858	EasyBlocker		


Rat EasyBlot® (GTX628474-01)



Goat EasyBlot® (GTX628547-01)



Sheep EasyBlot® (GTX628906-01)



International

Tel 886.3.6208988
Fax 886.3.6208989
Address 6F-2, No.89, Dongmei Rd., Hsinchu 300, Taiwan (R.O.C.)
E-mail international@genetex.com

USA

Toll-free 1.877.GeneTex(1.877.436.3839)
Fax 1.949.309.2888
Address 2456 Alton Parkway, Irvine, CA 92606 USA
E-mail info@genetex.com

