

Human Anti-Monkeypox Virus (MPXV) IgG ELISA Kit (H3L)

Overview

Catalog No. AK707018

Specifications H3L

Applications ELISA

Accession NP_536520.1

Stability and Storage When the kit was stored at the recommended temperature for 6 months, the signal intensity decreased by less than 20%. For unopened kits, if you want to prolong the storage time, please store the Standard, Detection A, Detection B and Microplate at - 20°C, the rest reagents should be store at 4°C.

Detection method Colorimetric

Sample type Plasma, Serum

Assay type Quantitative

Range 469 ng/mL - 30, 000 ng/mL

Recovery 80-120%

Shipping 2-8 °C

Note For Research Use Only.