

Human β Amyloid (1-40) ELISA KIT

Overview

Catalog No.	HY235028
Specifications	β Amyloid (1-40)
Applications	ELISA
Stability and Storage	The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 10% prior to the expiration date under appropriate storage condition.
Detection method	Colorimetric
Sample type	Plasma, Serum
Assay type	Quantitative
Sensitivity	200 pg/ml
Range	1-100ng/mL
Recovery	80-120%
Shipping	2-8 °C
Note	For research use only.