

# Imdevimab ELISA kit

## Overview

---

**Catalog No.** DK565038

**Specifications** Imdevimab

**Applications** ELISA

**Accession** P0DTC2

**Stability and Storage** When the kit was stored at the recommended temperature for 6 months, the signal intensity decreased by less than 20%. For unopened kits, if you want to prolong the storage time, please store the Standard, Detection A, Detection B and Microplate at - 20°C, the rest reagents should be store at 4°C.

**Detection method** Colorimetric

**Sample type** Plasma, Serum

**Assay type** Quantitative

**Range** 15.625 ng/mL - 1, 000 ng/mL

**Recovery** 80-120%

**Alternative name** REGN-10987, REGN10987, 2415933-40-1

**Shipping** 2-8 °C

**Note** For Research Use Only.