


Automatic ICE Maker Machine

UBT2101 Series



UBTECH™

Professional Creates the Quality, Innovation achieves the future

Snow Flake Icemaker- IMS Models



UBTECH Snow Flake Ice Maker machine IMS Models are compact and Durable machine with ice making capacity from 30kg/24h up to 500 kg/24h. The high quality stainless steel body ensures durability. It offers an increased efficiency up to 40% with energy saving up to 35%. Rapid flake shape ice production and inbuilt storage bin provides maximum performance under continuous operating condition.

Features

- Compact construction with high aesthetic appearance.
- Fully automatic operation system
- Air cooling, side ventilation system
- High efficiency economical CFC free compressor with low energy consumption system
- Over loading protection system, water shortage auto-detecting system

Application

UBTECH Snow Flake Ice Maker machines are widely used for biological, food preservation, biochemical and enzymatic studies in research and medical labs, health care facilities, food and beverage industries, pharmaceuticals, chemical industries, etc.



UNIVERSAL BIOTECHNOLOGY (NEW DELHI, INDIA)

Specification

Model No.	IMS-30	IMS-50	IMS-70	IMS-85	IMS-100
Ice production capacity	30kg/24h	50kg/24h	70kg/24h	85kg/24h	100kg/24h
Ice Storage Capacity	10kg	15kg	25kg	25kg	40kg
Power Supply	300W/220V	360W/220V	460W/220V	480W/220V	520W/220V
Outer Dimension	330×480×605mm	400×543×700mm	400×510×845mm	400×510×845mm	500×611×940mm
Water Inlet Mode	Tap water automatically inlet				
Cooling Mode	Air-cooling				
Refrigeration	R134a				
Ice Storage Capacity	Irregular small particles ice				

Model No.	IMS-150	IMS-200	IMS-250	IMS-300	IMS-500
Ice production capacity	150kg/24h	200kg/24h	250kg/24h	300kg/24h	500kg/24h
Ice Storage Capacity	40kg	60kg	60kg	200kg	300kg
Power Supply	680W/220V	1040W/220V	1100W/220V	1360W/220V	2100W/220V
Outer Dimension	500×611×940mm	680×740×1100mm	680×740×1100mm	560×870×1580mm	760×860×2050mm
Water Inlet Mode	Tap water automatically inlet				
Cooling Mode	Air-cooling				
Refrigeration	R134a				
Ice Storage Capacity	Irregular small particles ice				

*We can provide customize Model & Specification on special orders Request

UNIVERSAL BIOTECHNOLOGY (NEW DELHI, INDIA)