

Snow Icemaker- IMX Models



UBTECH Snow Ice Maker machine IMX Models are compact and Durable machine with ice making capacity from 60kg/24h up to 500 kg/24h. The high quality stainless steel body ensures durability. It offers an increased efficiency up to 40% with energy saving up to 35%. Rapid Snow ice production and inbuilt storage bin provides maximum performance under continuous operating condition.

Features

- Compact construction with high aesthetic appearance.
- Fully automatic operation system
- Air cooling, side ventilation system
- High efficiency economical CFC free compressor with low energy consumption system
- Over loading protection system, water shortage auto-detecting system

Application

UBTECH Snow Ice Maker machines are widely used for biological, food preservation, biochemical and enzymatic studies in research and medical labs, health care facilities, food and beverage industries, pharmaceuticals, chemical industries, etc.



IMX-100



IMX-500

Specification

Model No.	IMX-60	IMX-100	IMX-180	IMS-350	IMS-500
Ice production capacity	60kg/24h	100kg/24h	180kg/24h	350kg/24h	500kg/24h
Ice Storage Capacity	30kg	40kg	60kg	280kg	280kg
Power Supply	300W/220V	450W/220V	800W/220V	1200W/220V	1600W/220V
Outer Dimension	400×520×820mm	500×520×920mm	500×560×1220mm	760×830×1720mm	760×830×1720mm
Water Inlet Mode	Tap water automatically inlet				
Cooling Mode	Air-cooling / Water-cooling				
Refrigeration	R134a				
Ice Storage Capacity	Snow Ice				

UNIVERSAL BIOTECHNOLOGY (NEW DELHI, INDIA)