

# Inclacumab ELISA Kit

## Overview

---

**Catalog No.** DB885028

**Specifications** Inclacumab

**Applications** ELISA

**Accession** P16109

**Stability and Storage** When the kit was stored at the recommended temperature for 6 months, the signal intensity decreased by less than 20%. For unopened kits, if you want to prolong the storage time, please store the Standard, Detection A, Detection B and Microplate at - 20 °C, the rest reagents should be store at 4°C.

**Detection method** Colorimetric

**Sample type** Plasma, Serum

**Assay type** Quantitative

**Sensitivity** 5.67 ng/mL

**Range** 31.25 - 2, 000 ng/mL

**Recovery** 80-120%

**Alternative name** LC1004-002, RO4905417, 1256258-86-2

**Shipping** 2-8 °C

**Note** For Research Use Only.