

Ixekizumab ELISA Kit

Overview

Catalog No.	DS856038
Specifications	Ixekizumab
Applications	ELISA
Accession	Q16552
Stability and Storage	When the kit was stored at the recommended temperature for 6 months, the signal intensity decreased by less than 20%. For unopened kits, if you want to prolong the storage time, please store the Standard, Detection A, Detection B and Microplate at - 20°C, the rest reagents should be store at 4°C.
Detection method	Colorimetric
Sample type	Plasma, Serum
Assay type	Quantitative
Sensitivity	39.3 ng/mL
Range	78.13 - 5, 000 ng/mL
Recovery	80-120%
Alternative name	LY2439821, 1143503-69-8
Shipping	2-8 °C
Note	For Research Use Only.