

# Pipette Filler

## Levo Plus Pipette Filler

### Features

- Appropriate for most of the plastic and glass pipettes from 0.1 ~100ml
- Eight speeds of selection for aspiration and dispensing different liquids
- Large LCD display showing low battery warning and speed settings
- Enabling single-handed operation with minimum effort
- Light and ergonomic design providing easy usability
- High capacity Li-ion battery enable long runtime of operation
- Powerful pump fills a 25mL pipette in <5 seconds
- 0.45µm replaceable Hydrophobic filter
- Rechargeable during use

Large LCD display showing low battery warning and speed settings

Select the desired coarse speed by pressing the button +/-



Operating buttons fine-tuned for aspiration and dispensing

Operating button for motor and gravity dispensing



Hydrophobic filter 0.45µm

Silicone adapter for 0.1~100ml pipette tip, fully autoclavable



### Specifications

Aspirate Speeds	25ml < 5s (8 shift)
Dispense Speeds	Motor (8 shift) / Gravity
Battery	Lithium-Ion
Battery Service Life	More than 8 Hours of Intermittence Use
Charging Time	2~3 Hours
Pipette Types	Glass or Plastic Pipette (0.1~100mL), Pasteur Pipettes
Filter	0.45 µm Hydrophobic
Weight	200g

## TALK TO US

Toll Free: 1800 547 2314  
Office : +91 9910 73 8844  
Email : [info@unibiotech.in](mailto:info@unibiotech.in)