

# Magrolimab ELISA Kit

## Overview

---

**Catalog No.** DC359018

**Specifications** Magrolimab

**Applications** ELISA

**Accession** Q08722

**Stability and Storage** When the kit was stored at the recommended temperature for 6 months, the signal intensity decreased by less than 20%. For unopened kits, if you want to prolong the storage time, please store the Standard, Detection A, Detection B and Microplate at - 20 °C, the rest reagents should be store at 4°C.

**Detection method** Colorimetric

**Sample type** Plasma, Serum

**Assay type** Quantitative

**Sensitivity** 48.26 ng/mL

**Range** 62.5 - 4, 000 ng/mL

**Recovery** 80-120%

**Alternative name** Hu5F9-G4, 2169232-81-7

**Shipping** 2-8 °C

**Note** For Research Use Only.