

Monkeypox Virus (MPXV) A35R Antigen ELISA Kit

Overview

Catalog No. VK601018

Specifications A35R

Applications ELISA

Stability and Storage

When the kit was stored at the recommended temperature for 6 months, the signal intensity decreased by less than 20%. For unopened kits, if you want to prolong the storage time, please store the Standard, Detection A, Detection B and Microplate at - 20°C, the rest reagents should be store at 4°C.

Detection method Colorimetric

Sample type Plasma, Serum

Assay type Quantitative

Sensitivity 158.4 pg/mL

Range 375 pg/mL - 24, 000 pg/mL

Recovery 70-100%

Shipping 2-8 °C

Note For Research Use Only.