

Monkeypox Virus (MPXV) M1R Antigen ELISA Kit

Overview

Catalog No. VK785018

Specifications M1R

Applications ELISA

Stability and Storage

When the kit was stored at the recommended temperature for 6 months, the signal intensity decreased by less than 20%. For unopened kits, if you want to prolong the storage time, please store the Standard, Detection A, Detection B and Microplate at - 20°C, the rest reagents should be store at 4°C.

Detection method Colorimetric

Sample type Plasma, Serum

Assay type Quantitative

Sensitivity 1.9 pg/mL

Range 7.81 - 500 pg/mL

Recovery 85-120%

Shipping 2-8 °C

Note For Research Use Only.