

Monkeypox Virus (MPXV) A29L Antigen ELISA Kit-HS

Overview

Catalog No.	VK415028
Specifications	A29L
Applications	ELISA
Stability and Storage	The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 10% prior to the expiration date under appropriate storage condition.
Detection method	Colorimetric
Sample type	Plasma, Serum
Assay type	Quantitative
Sensitivity	3.93 pg/mL
Range	7.81 - 500 pg/mL
Recovery	80-120%
Shipping	2-8 °C
Note	For Research Use Only.