

# Nirsevimab ELISA Kit

## Overview

---

**Catalog No.** DK421018

**Specifications** Nirsevimab

**Applications** ELISA

**Accession** P03420

**Stability and Storage** When the kit was stored at the recommended temperature for 6 months, the signal intensity decreased by less than 20%. For unopened kits, if you want to prolong the storage time, please store the Standard, Detection A and Microplate at - 20 °C, the rest reagents should be store at 4°C.

**Detection method** Colorimetric

**Sample type** Plasma, Serum

**Assay type** Quantitative

**Sensitivity** 13.72 ng/mL

**Range** 46.88 - 3, 000 ng/mL

**Recovery** 80-120%

**Alternative name** MEDI-8897, MEDI8897, 1989556-22-0

**Shipping** 2-8 °C

**Note** For Research Use Only.