

# Olaratumab ELISA Kit

## Overview

---

<b>Catalog No.</b>	DB927018
<b>Specifications</b>	Olaratumab
<b>Applications</b>	ELISA
<b>Accession</b>	P16234
<b>Stability and Storage</b>	When the kit was stored at the recommended temperature for 6 months, the signal intensity decreased by less than 20%. For unopened kits, if you want to prolong the storage time, please store the Standard, Detection A, Detection B and Microplate at - 20°C, the rest reagents should be store at 4°C.
<b>Detection method</b>	Colorimetric
<b>Sample type</b>	Plasma, Serum
<b>Assay type</b>	Quantitative
<b>Sensitivity</b>	112.17 ng/mL
<b>Range</b>	156.25 - 10, 000 ng/mL
<b>Recovery</b>	80-120%
<b>Alternative name</b>	3G3, IMC-3G3, LY3012207, 1024603-93-7
<b>Shipping</b>	2-8 °C
<b>Note</b>	For Research Use Only.