

Pembrolizumab ELISA Kit

Overview

Catalog No.	DS870048
Specifications	Pembrolizumab
Applications	ELISA
Accession	Q15116
Stability and Storage	The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 10% prior to the expiration date under appropriate storage condition.
Detection method	Colorimetric
Sample type	Plasma, Serum
Assay type	Quantitative
Sensitivity	0.156 µg/ml
Range	0.31-5 µg/mL
Recovery	80-120%
Alternative name	MK-3475
Shipping	2-8 °C
Note	For Research Use Only.