

TNM-FH medium for insect Cells

Cat. No: PM152010

Size: 500mL

General Information

| | |
|--------------------------|---------------------|
| Product Form | Liquid |
| Concentration | 1× |
| pH | 6.1-6.4 |
| D-Glucose | 700mg/L |
| HEPES | Negative |
| L-Glutamine | 600mg/L |
| NaHCO₃ | 350mg/L |
| Phenol red | Negative |
| Storage | 2-8℃ ,Shading Light |
| Shipping | RT |
| Expiration date | 24 months |

Background

TNM-FH medium is an improved Grace's insect medium. After proper supplementation, it can be used for the culture of various Lepidoptera insect cells.

TNM-FH medium contains amino acids, vitamins, inorganic salts and other ingredients required for various types of cell culture, but does not contain protein, lipids or any growth factors, so this product needs to be used with serum or serum-free additives.

Notes

1. This product is for research use only.;
2. This product is sterilized by 0.1μm filtration.;
3. It is necessary to pay attention to the aseptic operation and avoid the contamination during the culture.