

Romiplostim ELISA Kit

Overview

Catalog No. YB762018

Specifications Romiplostim

Applications ELISA

Accession CAS: 267639-76-9

Stability and Storage When the kit was stored at the recommended temperature for 12 months, the signal intensity decreased by less than 20%. For unopened kits, if you want to prolong the storage time, please store the Standard, Detection A, Detection B and Microplate at - 20 °C, the rest reagents should be store at 4°C.

Detection method Colorimetric

Sample type Plasma, Serum

Assay type Quantitative

Sensitivity 11.31 pg/mL

Range 11.72 - 750 pg/mL

Recovery 80-120%

Alternative name 267639-76-9, romiplostim (Nplate)

Shipping 2-8 °C

Note For Research Use Only.