

# Romosozumab ELISA Kit

## Overview

---

**Catalog No.** DV185028

**Specifications** Romosozumab

**Applications** ELISA

**Accession** Q9BQB4

**Stability and Storage** When the kit was stored at the recommended temperature for 6 months, the signal intensity decreased by less than 20%. For unopened kits, if you want to prolong the storage time, please store the Standard, Detection A, Detection B and Microplate at - 20°C, the rest reagents should be store at 4°C.

**Detection method** Colorimetric

**Sample type** Plasma, Serum

**Assay type** Quantitative

**Sensitivity** 42.86 ng/mL

**Range** 125 - 8, 000 ng/mL

**Recovery** 80-120%

**Alternative name** 785A070802, AMG 785, CDP-7851, romosozumab-aqqg, sclerostin Ab, 909395-70-6

**Shipping** 2-8 °C

**Note** For Research Use Only.