

Rosopatamab ELISA Kit

Overview

Catalog No. DC669028

Specifications Rosopatamab

Applications ELISA

Accession Q04609

Stability and Storage When the kit was stored at the recommended temperature for 6 months, the signal intensity decreased by less than 20%. For unopened kits, if you want to prolong the storage time, please store the Standard, Detection A, Detection B and Microplate at - 20 °C, the rest reagents should be store at 4°C.

Detection method Colorimetric

Sample type Plasma, Serum

Assay type Quantitative

Sensitivity 63.23 ng/mL

Range 468.75 - 30, 000 ng/mL

Recovery 80-120%

Alternative name J-591, MLN-591, huJ-591, Huj591-DOTA, Rosopatamab tetraxetan, 2260767-49-3

Shipping 2-8 °C

Note For Research Use Only.