

Sacituzumab ELISA Kit

Overview

Catalog No.	DY373018
Specifications	Sacituzumab
Applications	ELISA
Accession	P09758
Stability and Storage	When the kit was stored at the recommended temperature for 6 months, the signal intensity decreased by less than 20%. For unopened kits, if you want to prolong the storage time, please store the Standard, Detection A, Detection B and Microplate at - 20°C, the rest reagents should be store at 4°C.
Detection method	Colorimetric
Sample type	Plasma, Serum
Assay type	Quantitative
Sensitivity	40.3 ng/mL
Range	78.13 - 5, 000 ng/mL
Recovery	80-120%
Alternative name	IMMU-132, hRS7-SN-38, hRS7-SN-38-ADC, hRS7-[CL-SN-38], , 1796566-95-4
Shipping	2-8 °C
Note	For Research Use Only.