

SARS-CoV-2 Surrogate Virus Neutralization Test Kit (JN.1)

Overview

| | |
|------------------------------|---|
| Catalog No. | AK565098 |
| Specifications | SARS-CoV-2 Neutralization Antibody |
| Applications | ELISA |
| Stability and Storage | The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 10% prior to the expiration date under appropriate storage condition. |
| Detection method | Colorimetric |
| Sample type | Plasma, Serum |
| Assay type | Quantitative |
| Sensitivity | 16.81 ng/mL |
| Range | 31.25 - 2000 ng/mL |
| Recovery | 80-120% |
| Shipping | 2-8 °C |
| Note | For research use only. |