

Sotrovimab ELISA kit

Overview

Catalog No.	DK565018
Specifications	Sotrovimab
Applications	ELISA
Accession	P0DTC2
Stability and Storage	When the kit was stored at the recommended temperature for 6 months, the signal intensity decreased by less than 20%. For unopened kits, if you want to prolong the storage time, please store the Standard, Detection A, Detection B and Microplate at - 20°C, the rest reagents should be store at 4°C.
Detection method	Colorimetric
Sample type	Plasma, Serum
Assay type	Quantitative
Range	156.25 ng/mL - 10, 000 ng/mL
Recovery	80-120%
Alternative name	VIR-7831, GSK4182136, S309, GSK-4182136, 2423014-07-5
Shipping	2-8 °C
Note	For Research Use Only.