

Urelumab ELISA Kit

Overview

Catalog No. DC336018

Specifications Urelumab

Applications ELISA

Accession Q07011

Stability and Storage When the kit was stored at the recommended temperature for 6 months, the signal intensity decreased by less than 20%. For unopened kits, if you want to prolong the storage time, please store the Standard, Detection A, Detection B and Microplate at - 20°C, the rest reagents should be store at 4°C.

Detection method Colorimetric

Sample type Plasma, Serum

Assay type Quantitative

Sensitivity 18.73 ng/mL

Range 31.25 - 2, 000 ng/mL

Recovery 80-120%

Alternative name BMS-663513, 934823-49-1

Shipping 2-8 °C

Note For Research Use Only.