





Cytokines

ELISA Kits for Cytokines

• Cytokines

Cytokine is a class of small proteins (5-20kDa) that are important in cell signaling. Their release has an effect on the behaviour of cells around them. Cytokines are involved in autocrine signaling, paracrine signaling and endocrine signaling as immunomodulating agents. Their definite distinction from hormones is still an ongoing research. Cytokines include chemokines, interferons, interleukins, lymphokines and tumor necrosis factors but generally not hormones or growth factors (despite some overlap in the terminology). Cytokines are produced by a broad range of cells, including immune cells like macrophages, B lymphocytes, T lymphocytes and mast cells, as well as endothelial cells, fibroblasts and various stromal cells.

Cytokine is important in health and disease research, specifically in host response to infection, immune response, inflammation, trauma, sepsis, cancer and reproduction.



Chemokines

Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-H5490	Human CXCR4(Chemokine C-X-C-Motif Receptor) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich
E-EL-H0679	Human CCR1 (CC-Chemokine Receptor 1) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-H5489	Human CCR4(Chemokine C-C-Motif Receptor 4) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich
E-EL-H0023	Human RANTES (Regulated On Activation, Normal T-Cell Expressed and Secreted) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-H0052	Human SDF-1 (Stromal Cell Derived Factor 1) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-H0854	Human CXCR3 (Chemokine CXC Receptor 3) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-H1877	Human CCR2 (Chemokine C-C-Motif Receptor 2) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-H0020	Human MCP-1 (Monocyte Chemotactic Protein 1) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-H0704	Human CCR8 (Chemokine C-C-Motif Receptor 8) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-H1528	Human MIP-3 β /ELC/CCL19 (Macrophage Inflammatory Protein 3 β) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-H0022	Human MIP-1 β (Macrophage Inflammatory Protein 1 Beta) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-H0026	Human TARC (Thymus Activation Regulated Chemokine) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-H0028	Human SLC (Secondary Lymphoid Tissue Chemokine) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich
E-EL-H0050	Human IP-10/CXCL10 (Interferon Gamma Induced Protein 10kDa) ELISA Kit	7.81--500pg/mL	4.69pg/mL	Sandwich
E-EL-H0029	Human MDC (Macrophage Derived Chemokine) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-H0048	Human IL-8 (Interleukin 8) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-H0707	Human CX3CR1 (Chemokine C-X3-C-Motif Receptor 1) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-H2578	Human CXCR2(CXC-Chemokine Receptor 2) ELISA Kit	0.78--50ng/mL	0.47ng/mL	Sandwich
E-EL-H0027	Human MIP-3 α (Macrophage Inflammatory Protein 3 Alpha) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-H0025	Human ECF/CCL11 (Eosinophil Chemotactic Factor) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-H0021	Human MIP-1 α (Macrophage Inflammatory Protein 1 Alpha) ELISA Kit	23.44--1500pg/mL	14.06pg/mL	Sandwich
E-EL-H0053	Human BLC (B-Lymphocyte Chemoattractant) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-H1009	Human Mly/CXCL9 (Monocyte Interferon Gamma Inducing Factor) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-H0046	Human ENA-78 (Epithelial Neutrophil Activating Peptide 78) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-H2144	Human TECK/CCL25 (Thymus Expressed Chemokine) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-H1904	Human GRO β (Growth Regulated Oncogene Beta) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-H0051	Human I-TAC (Interferon Inducible T-cell Alpha Chemoattractant) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-H0699	Human CCL1 (Chemokine C-C-Motif Ligand 1) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich

Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-H1691	Human Lptn/XCL1 (Lymphotactin) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-H1159	Human MCP-4/CCL13 (Monocyte Chemotactic Protein 4) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-H1158	Human MCP-2/CCL8 (Monocyte Chemotactic Protein 2) ELISA Kit	6.25--400pg/mL	3.75pg/mL	Sandwich
E-EL-H0047	Human GCP-2 (Granulocyte Chemotactic Protein 2) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-H0032	Human CTACK (Cutaneous T-cell Attracting Chemokine) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-H0962	Human XCR1 (XC-Chemokine Receptor 1) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-H0045	Human GRO α /CXCL1 (Growth Regulated Oncogene Alpha) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-H1270	Human PARC/CCL18 (Pulmonary Activation Regulated Chemokine) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich
E-EL-H0054	Human CXCL16 (Chemokine C-X-C-Motif Ligand 16) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-H0853	Human CXCR1 (CXC-Chemokine Receptor 1) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-H0031	Human MIP4 α (Macrophage Inflammatory Protein 4 Alpha) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-H1348	Human MPIF-1/CCL23 (Myeloid Progenitor Inhibitory Factor 1) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-H2282	Human MEC/CCL28 (Mucosae Associated Epithelia Chemokine) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-H0030	Human MPIF2 (Myeloid Progenitor Inhibitory Factor 2) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich
E-EL-H0044	Human CX3CL1 (Chemokine C-X3-C-Motif Ligand 1) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-H0701	Human CCL16 (Chemokine C-C-Motif Ligand 16) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-H0600	Human BRAK/CXCL14 (Breast and Kidney Expressed Chemokine) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-H0700	Human CCL14 (Chemokine C-C-Motif Ligand 14) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-H0702	Human CCL3L1 (Chemokine C-C-Motif Ligand 3 Like Protein 1) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich
E-EL-H1905	Human GRO γ (Growth Regulated Oncogene Gamma) ELISA Kit	7.81--500pg/mL	4.69pg/mL	Sandwich
E-EL-H0703	Human CCL4L1 (Chemokine C-C-Motif Ligand 4 Like Protein 1) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-H0714	Human CKLF (Chemokine Like Factor) ELISA Kit	125--8000pg/mL	75pg/mL	Sandwich
E-EL-H2216	Human PF4 (Platelet Factor 4) ELISA Kit	1.56--100ng/mL	0.94ng/mL	Sandwich
E-EL-M0009	Mouse RANTES (Regulated On Activation, Normal T-Cell Expressed and Secreted) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-M0373	Mouse CXCR3 (CXC-chemokine receptor 3) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-M0006	Mouse MCP-1 (Monocyte Chemotactic Protein 1) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M1094	Mouse SDF-1/CXCL12(Stromal Cell Derived Factor 1) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-M0008	Mouse MIP-1 β (Macrophage Inflammatory Protein 1 Beta) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0265	Mouse CCR8(Chemokine C-C-Motif Receptor 8) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-M0021	Mouse IP-10/CXCL10 (Interferon Gamma Induced Protein 10kDa) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich

Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-M0007	Mouse MIP-1 α (Macrophage Inflammatory Protein 1 Alpha) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0769	Mouse MIP-3 β (Macrophage Inflammatory Protein 3 β) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0012	Mouse TARC (Thymus Activation Regulated Chemokine) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-M0145	Mouse SLC (Secondary Lymphoid Tissue Chemokine) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0014	Mouse MDC/CCL22 (Macrophage-Derived Chemokine) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0013	Mouse MIP-3 α (Macrophage Inflammatory Protein 3 Alpha) ELISA Kit	7.81--500pg/mL	4.69pg/mL	Sandwich
E-EL-M0464	Mouse ECF/CCL11 (Eosinophil Chemotactic Factor) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-M0268	Mouse CX3CR1 (Chemokine C-X3-C-Motif Receptor 1) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-M0020	Mouse CXCL9/MIG(Monokine induced by interferon-gamma) ELISA Kit	46.88--3000pg/mL	28.16pg/mL	Sandwich
E-EL-M0056	Mouse I-TAC (Interferon Inducible T-Cell Alpha Chemoattractant) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-M0018	Mouse GRO α /CXCL1 (Growth Regulated Oncogene Alpha) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0190	Mouse BLC (B-Lymphocyte Chemoattractant) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0471	Mouse ENA-78 (Epithelial Neutrophil Activating Peptide 78) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M1148	Mouse TECK (Thymus Expressed Chemokine) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-M0263	Mouse CCL1 (Chemokine C-C-Motif Ligand 1) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-M0794	Mouse MCP-2 (Monocyte Chemotactic Protein 2) ELISA Kit	6.25--400pg/mL	3.75pg/mL	Sandwich
E-EL-M0796	Mouse MCP-4 (Monocyte Chemotactic Protein 4) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-M0991	Mouse PARC (Pulmonary Activation Regulated Chemokine) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich
E-EL-M0766	Mouse MIP-1 δ (Macrophage Inflammatory Protein 1 δ) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-M0763	Mouse Lptn/XCL1 (Lymphotactin) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-M0267	Mouse CX3CL1 (Chemokine C-X3-C-Motif Ligand 1) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-M0371	Mouse CTACK (Cutaneous T-Cell Attracting Chemokine) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-M0466	Mouse MPIF2 (Myeloid Progenitor Inhibitory Factor 2) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0592	Mouse GCP-2 (Granulocyte Chemotactic Protein 2) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-M0270	Mouse CXCL16 (Chemokine C-X-C-Motif Ligand 16) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-M0147	Mouse GRO γ (Growth Regulated Oncogene Gamma) ELISA Kit	7.81--500pg/mL	4.69pg/mL	Sandwich
E-EL-M0206	Mouse BRAK (Breast And Kidney Expressed Chemokine) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-M0264	Mouse CCL3L1 (Chemokine C-C-Motif Ligand 3 Like Protein 1) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich
E-EL-M0926	Mouse PF4 (Platelet Factor 4) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich
E-EL-M0271	Mouse CKLF (Chemokine Like Factor) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich

Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-M0048	Mouse IFN- γ (Interferon Gamma) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0049	Mouse TNF- α (Tumor Necrosis Factor Alpha) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-R0845	Rat RANTES (Regulated On Activation, Normal T-Cell Expressed and Secreted) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0287	Rat CXCR3 (CXC-Chemokine Receptor 3) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-R0633	Rat MCP-1 (Monocyte Chemotactic Protein 1) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-R2581	Rat CCR1 (CC-Chemokine Receptor 1) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-R0922	Rat SDF-1 (Stromal Cell Derived Factor 1) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-R0603	Rat MIP-1 β (Macrophage Inflammatory Protein 1 Beta) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0185	Rat CCR8 (Chemokine C-C-Motif Receptor 8) ELISA Kit	1.56--100ng/mL	0.94ng/mL	Sandwich
E-EL-R0546	Rat IP-10/CXCL10 (Interferon Gamma Induced Protein 10kDa) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-R0602	Rat MIP-1 α (Macrophage Inflammatory Protein 1 Alpha) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-R0606	Rat MIP-3 β /ELC/CCL19(Macrophage Inflammatory Protein 3 β) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0969	Rat TARC (Thymus Activation Regulated Chemokine) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-R0627	Rat SLC (Secondary Lymphoid Tissue Chemokine) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-R0609	Rat MDC/CCL22 (Macrophage-Derived Chemokine) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0605	Rat MIP-3 α (Macrophage Inflammatory Protein 3 Alpha) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0368	Rat ECF/CCL11(Eosinophil Chemotactic Factor) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-R0188	Rat CX3CR1 (Chemokine C-X3-C-Motif Receptor 1) ELISA Kit	1.56--100ng/mL	0.94ng/mL	Sandwich
E-EL-R0854	Rat Mly/CXCL9 (Monocyte Interferon Gamma Inducing Factor) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-R2433	Rat I-TAC (Interferon Inducible T-Cell Alpha Chemoattractant) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-R0696	Rat GRO β (Growth Regulated Oncogene Beta) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-R0003	Rat GRO α /CXCL1 (Growth Regulated Oncogene Alpha) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0112	Rat BLC (B-Lymphocyte Chemoattractant) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0370	Rat ENA-78 (Epithelial Neutrophil Activating Peptide 78) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0970	Rat TECK (Thymus Expressed Chemokine) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-R0184	Rat CCL1 (Chemokine C-C-Motif Ligand 1) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-R0600	Rat Lptn/XCL1 (Lymphotactin) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-R0385	Rat CX3CL1 (Chemokine C-X3-C-Motif Ligand 1) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-R0285	Rat CTACK (Cutaneous T-Cell Attracting Chemokine) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich

Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-R0356	Rat MPIF2 (Myeloid Progenitor Inhibitory Factor 2) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-R0453	Rat GCP-2 (Granulocyte Chemotactic Protein 2) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-R0189	Rat CXCL16 (Chemokine C-X-C-Motif Ligand 16) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-R0697	Rat GROγ (Growth Regulated Oncogene Gamma) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-R0128	Rat BRAK (Breast And Kidney Expressed Chemokine) ELISA Kit	93.75--6000pg/mL	56.25pg/mL	Sandwich
E-EL-R0190	Rat CKLF (Chemokine Like Factor) ELISA Kit	7.81--500pg/mL	4.69pg/mL	Sandwich
E-EL-R0009	Rat IFN-γ (Interferon Gamma) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-R0019	Rat TNF-α (Tumor Necrosis Factor Alpha) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich

Interferons

Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-H0108	Human IFN-γ (Interferon Gamma) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-H0085	Human IFN-β (Interferon Beta) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-H2491	Human IRF3 (Interferon Regulatory Factor 3) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-H2130	Human STAT1 (Signal Transducer and Activator of Transcription 1) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-H1545	Human IRF5 (Interferon Regulatory Factor 5) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-H2131	Human STAT2 (Signal Transducer and Activator of Transcription 2) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-H0102	Human IL-6 (Interleukin 6) ELISA Kit	7.81--500pg/mL	4.69pg/mL	Sandwich
E-EL-H0253	Human IL-18 (Interleukin 18) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-H0203	Human ISG15 (Interferon Stimulated Gene 15) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-H0099	Human IL-2 (Interleukin 2) ELISA Kit	125--8000pg/mL	75pg/mL	Sandwich
E-EL-H5552	Human JAK1 (Janus Kinase 1) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich

Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-H1653	Human Tyk-2 (Tyrosine Kinase 2) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-H0103	Human IL-10 (Interleukin 10) ELISA Kit	7.81--500pg/mL	4.69pg/mL	Sandwich
E-EL-M0713	Mouse IFN- α (Interferon Alpha) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0048	Mouse IFN- γ (Interferon Gamma) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0033	Mouse IFN- β (Interferon Beta)ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0969	Mouse PIAS1 (Protein Inhibitor Of Activated STAT 1) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-M2449	Mouse IRF5 (Interferon Regulatory Factor 5)	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-M2688	Mouse IFN α 4(Interferon Alpha 4) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-M1056	Mouse STAT2 (Signal Transducer And Activator Of Transcription 2) ELISA Kit	0.63--40ng/mL	0.38ng/mL	Sandwich
E-EL-M0730	Mouse IL-18 (Interleukin 18) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0044	Mouse IL-6 (Interleukin 6) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-M0049	Mouse TNF- α (Tumor Necrosis Factor Alpha) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-M0042	Mouse IL-2 (Interleukin 2) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0046	Mouse IL-10 (Interleukin 10) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0715	Mouse IFI16/p16 (Interferon-inducible Protein 16) ELISA Kit	23.44--1500pg/mL	14.06pg/mL	Sandwich
E-EL-M0021	Mouse IP-10/CXCL10 (Interferon Gamma Induced Protein 10kDa) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-R0009	Rat IFN- γ (Interferon Gamma) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-R0545	Rat IFN- β (Interferon Beta) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R1131	Rat IRF4(Interferon Regulatory Factor 4) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-R0567	Rat IL-18 (Interleukin 18) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0015	Rat IL-6 (Interleukin 6) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-R0019	Rat TNF- α (Tumor Necrosis Factor Alpha) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich
E-EL-R0013	Rat IL-2 (Interleukin 2) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-R0016	Rat IL-10 (Interleukin 10) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-R0547	Rat IFI16/p16 (Interferon-Inducible Protein 16) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R1227	Rat JAK2 (Janus Kinase 2) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-R0546	Rat IP-10/CXCL10 (Interferon Gamma Induced Protein 10kDa) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-H0102	Human IL-6 (Interleukin 6) ELISA Kit	7.81--500pg/mL	4.69pg/mL	Sandwich
E-EL-H0149	Human IL-1 β (Interleukin 1 Beta) ELISA Kit	7.81--500pg/mL	4.69pg/mL	Sandwich
E-EL-H0089	Human IL-1ra/IL-1F3 (Interleukin 1 Receptor Antagonist) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich

Interleukins

Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-H0103	Human IL-10 (Interleukin 10) ELISA Kit	7.81--500pg/mL	4.69pg/mL	Sandwich
E-EL-H1030	Human IL-2R α (Interleukin 2 Receptor Alpha) ELISA Kit	93.75--6000pg/mL	56.25pg/mL	Sandwich
E-EL-H0099	Human IL-2 (Interleukin 2) ELISA Kit	125--8000pg/mL	75pg/mL	Sandwich
E-EL-H1028	Human IL1R1 (Interleukin 1 Receptor Type I) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-H0088	Human IL-1 α (Interleukin 1 Alpha) ELISA Kit	3.91--250pg/mL	2.34pg/mL	Sandwich
E-EL-H0048	Human IL-8 (Interleukin 8) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-H0101	Human IL-4 (Interleukin 4) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-H0253	Human IL-18 (Interleukin 18) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-H1031	Human IL-2R β (Interleukin 2 Receptor Beta) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich
E-EL-H0191	Human IL-5 (Interleukin 5) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-H0192	Human IL-6R (Interleukin 6 Receptor) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-H0090	Human IL-1R2 (Interleukin 1 Receptor Type II) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-H0104	Human IL-13(Interleukin 13)ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-H0082	Human sgp130 (soluble Glucoprotein 130) ELISA Kit	0.256--16ng/mL	0.144ng/mL	Sandwich
E-EL-H0100	Human IL-3 (Interleukin 3) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-H5022	Human IL-11 (Interleukin 11) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-H1615	Human sST2 (soluble ST2)ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-H0648	Human IL-7 (Interleukin 7) ELISA Kit	7.81--500pg/mL	4.69pg/mL	Sandwich
E-EL-H0109	Human TNF- α (Tumor Necrosis Factor Alpha) ELISA Kit	7.81--500pg/mL	4.69pg/mL	Sandwich
E-EL-H0222	Human IL-15 (Interleukin 15) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-H0180	Human IL-9 (Interleukin 9) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-H1647	Human IL-12 p35 (Interleukin 12 p35) ELISA Kit	7.81--500pg/mL	4.69pg/mL	Sandwich
E-EL-H0235	Human IL-16 (Interleukin 16) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-H2338	Human IL-27 (Interleukin 27) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-H0424	Human IL10R β (Interleukin 10 Receptor Beta) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-H2578	Human CXCR2(CXC-Chemokine Receptor 2) ELISA Kit	0.78--50ng/mL	0.47ng/mL	Sandwich
E-EL-H0106	Human IL-22 (Interleukin 22) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-H0016	Human CASP1 (Caspase 1) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich

Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-H2135	Human STAT6 (Signal Transducer and Activator of Transcription 6) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-H2571	Human IL-37 (Interleukin 37) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-H1648	Human IL-25 (Interleukin 25) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-H1662	Human IL18BP (Interleukin 18 Binding Protein) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-H1002	Human IRAK3 (Interleukin 1 Receptor Associated Kinase 3) ELISA Kit	0.78--50ng/mL	0.47ng/mL	Sandwich
E-EL-H0254	Human IL-19 (Interleukin 19) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-H0279	Human IL-20 (Interleukin 20) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-H2450	Human IL-21 (Interleukin 21) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-H0108	Human IFN- γ (Interferon Gamma) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-H0853	Human CXCR1 (CXC-Chemokine Receptor 1) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-H0216	Human IL-32(Interleukin 32)ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-H2132	Human STAT3 (Signal Transducer and Activator of Transcription 3) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-H1650	Human IL-34 (Interleukin 34) ELISA Kit	93.75--6000pg/mL	56.25pg/mL	Sandwich
E-EL-H0081	Human GM-CSF (Granulocyte Macrophage Colony Stimulating Factor) ELISA Kit	7.81--500pg/mL	4.69pg/mL	Sandwich
E-EL-H0043	Human CRP (C-Reactive Protein) ELISA Kit	0.39--25ng/mL	0.23ng/mL	Sandwich
E-EL-H0361	Human CSF2R β (Colony Stimulating Factor 2 Receptor Beta) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-H2247	Human OSM (Oncostatin M) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-H5559	Human JAK3 (Janus Kinase 3) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich
E-EL-H1032	Human NFIL3 (Nuclear Factor, Interleukin 3 Regulated) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-M0044	Mouse IL-6 (Interleukin 6) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-M0037	Mouse IL-1 β (Interleukin 1 Beta) ELISA Kit	7.81--500pg/mL	4.69pg/mL	Sandwich
E-EL-M0038	Mouse IL-1ra (Interleukin 1 Receptor Antagonist) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-M0046	Mouse IL-10 (Interleukin 10) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0042	Mouse IL-2 (Interleukin 2) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0036	Mouse IL-1 α (Interleukin 1 Alpha)ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0043	Mouse IL-4 (Interleukin 4) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-M0730	Mouse IL-18 (Interleukin 18) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0722	Mouse IL-5 (Interleukin 5) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M2453	Mouse IL-6R (Interleukin 6 Receptor) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-M2451	Mouse IL-12/p40 (Interleukin 12 p40) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich

Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-M0727	Mouse IL-13 (Interleukin 13) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-M0721	Mouse IL-3 (Interleukin 3) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0049	Mouse TNF- α (Tumor Necrosis Factor Alpha) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-M0560	Mouse gp130 (Glucoprotein 130) ELISA Kit	125--8000pg/mL	75pg/mL	Sandwich
E-EL-M2443	Mouse IL-1R2 (Interleukin 1 Receptor Type II) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich
E-EL-M0047	Mouse IL-17 (Interleukin 17) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0723	Mouse IL-7 (Interleukin 7) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0728	Mouse IL-15 (Interleukin 15) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M2450	Mouse IL-12 p35 (Interleukin 12 p35) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0724	Mouse IL-9 (Interleukin 9) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0048	Mouse IFN- γ (Interferon Gamma) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0729	Mouse IL-16 (Interleukin 16) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0732	Mouse IL-27 (Interleukin 27) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-M2708	Mouse IL17RA(Interleukin-17 receptor A) ELISA kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M2446	Mouse IL-22 (Interleukin 22) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0201	Mouse CASP1 (Caspase 1) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-M1273	Mouse IL-21 (Interleukin 21) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-M1061	Mouse STAT6 (Signal Transducer And Activator Of Transcription 6) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich
E-EL-M0032	Mouse GM-CSF (Granulocyte-Macrophage Colony Stimulating Factor) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0739	Mouse IL18BP (Interleukin 18 Binding Protein) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-M2455	Mouse IL-20 (Interleukin 20) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-M0187	Mouse IL-25 (Interleukin 25) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0053	Mouse CRP (C-Reactive Protein) ELISA Kit	0.39--25ng/mL	0.23ng/mL	Sandwich
E-EL-M0051	Mouse TGF- β 1 (Transforming Growth Factor β 1) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-M0737	Mouse IL-32 (Interleukin 32) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0040	Mouse LIF (Leukemia Inhibitory Factor) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0859	Mouse OSM (Oncostatin-M) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0006	Mouse MCP-1 (Monocyte Chemotactic Protein 1) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0720	Mouse IL-34 (Interleukin 34) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-M0636	Mouse SCF (Stem Cell Factor) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich

Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-M0077	Mouse sICAM-1/CD54 (Soluble Intercellular Adhesion Molecule 1) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-M2687	Mouse MyD88 (Myeloid Differentiation Factor 88) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich
E-EL-M0031	Mouse G-CSF (Granulocyte Colony Stimulating Factor 3) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0713	Mouse IFN- α (Interferon Alpha) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0837	Mouse NFIL3 (Nuclear Factor, Interleukin 3 Regulated) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-R0015	Rat IL-6 (Interleukin 6) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-R0012	Rat IL-1 β (Interleukin 1 Beta) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-R0552	Rat IL-1ra (Interleukin 1 Receptor Antagonist) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich
E-EL-R0016	Rat IL-10 (Interleukin 10) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-R0013	Rat IL-2 (Interleukin 2) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-R0011	Rat IL-1 α (Interleukin 1 Alpha) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0014	Rat IL-4 (Interleukin 4) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0567	Rat IL-18 (Interleukin 18) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0558	Rat IL-5 (Interleukin 5) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0896	Rat IL-6R (Interleukin 6 Receptor) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich
E-EL-R0063	Rat IL-12/p40 (Interleukin 12 p40) ELISA Kit	12.5--800pg/mL	7.5pg/mL	Sandwich
E-EL-R0563	Rat IL-13 (Interleukin 13) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0556	Rat IL-3 (Interleukin 3) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0019	Rat TNF- α (Tumor Necrosis Factor Alpha) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich
E-EL-R0427	Rat gp130 (Glucoprotein 130) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-R2437	Rat IL-1R2 (Interleukin 1 Receptor Type II) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich
E-EL-R0566	Rat IL-17 (Interleukin 17) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0559	Rat IL-7 (Interleukin 7) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0564	Rat IL-15 (Interleukin 15) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R2442	Rat IL-12 p35 (Interleukin 12 p35) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0561	Rat IL-9 (Interleukin 9) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0009	Rat IFN- γ (Interferon Gamma) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-R0565	Rat IL-16 (Interleukin 16) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R1223	Rat IL-27 (Interleukin 27) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R2440	Rat IL-22 (Interleukin 22) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich

Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-R0371	Rat CASP1 (Caspase 1) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-R0568	Rat IL-21 (Interleukin 21) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0008	Rat GM-CSF (Granulocyte-Macrophage Colony Stimulating Factor) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R2426	Rat IL18BP (Interleukin 18 Binding Protein) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-R2446	Rat IL-20 (Interleukin 20) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-R2394	Rat IL-25 (Interleukin 25) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0022	Rat CRP (C-Reactive Protein) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-R0084	Rat TGF- β 1 (Transforming Growth Factor β 1) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-R0571	Rat IL-32 (Interleukin 32) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0585	Rat LIF (Leukemia Inhibitory Factor) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0692	Rat OSM (Oncostatin-M) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0633	Rat MCP-1 (Monocyte Chemotactic Protein 1) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-R2423	Rat IL-34 (Interleukin 34) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-R0919	Rat SCF (Stem Cell Factor) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R1227	Rat JAK2 (Janus Kinase 2) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-R0454	Rat G-CSF (Granulocyte Colony Stimulating Factor) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0672	Rat NFIL3 (Nuclear Factor, Interleukin 3 Regulated) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich

Lymphokines

Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-H0099	Human IL-2 (Interleukin 2) ELISA Kit	125--8000pg/mL	75pg/mL	Sandwich
E-EL-H0101	Human IL-4 (Interleukin 4) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-H0108	Human IFN- γ (Interferon Gamma) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-H1817	Human GPI (Glucose 6 Phosphate Isomerase) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-H1530	Human MIF (Macrophage Migration Inhibitory Factor) ELISA Kit	0.78--50ng/mL	0.47ng/mL	Sandwich
E-EL-H1030	Human IL-2R α (Interleukin 2 Receptor Alpha) ELISA Kit	93.75--6000pg/mL	56.25pg/mL	Sandwich
E-EL-H0191	Human IL-5 (Interleukin 5) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-H1031	Human IL-2R β (Interleukin 2 Receptor Beta) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich
E-EL-M0042	Mouse IL-2 (Interleukin 2) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0043	Mouse IL-4 (Interleukin 4) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-M0048	Mouse IFN- γ (Interferon Gamma) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0722	Mouse IL-5 (Interleukin 5) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0565	Mouse GPI (Glucose 6 Phosphate Isomerase) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-M0721	Mouse IL-3 (Interleukin 3) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0013	Rat IL-2 (Interleukin 2) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-R0014	Rat IL-4 (Interleukin 4) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0009	Rat IFN- γ (Interferon Gamma) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-R0558	Rat IL-5 (Interleukin 5) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0556	Rat IL-3 (Interleukin 3) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich

Tumor necrosis factors

Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-H0109	Human TNF- α (Tumor Necrosis Factor Alpha) ELISA Kit	7.81--500pg/mL	4.69pg/mL	Sandwich
E-EL-H0217	Human TNFRSF1A (Tumor Necrosis Factor Receptor Superfamily, Member 1A) ELISA Kit	7.81--500pg/mL	4.69pg/mL	Sandwich
E-EL-H2436	Human TNFRSF1B (Tumor Necrosis Factor Receptor Superfamily, Member 1B) ELISA Kit	7.81--500pg/mL	4.69pg/mL	Sandwich
E-EL-H0102	Human IL-6 (Interleukin 6) ELISA Kit	7.81--500pg/mL	4.69pg/mL	Sandwich
E-EL-H0111	Human VEGF-A (Vascular Endothelial Cell Growth Factor A) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-H0086	Human IGF-1 (Insulin-Like Growth Factor 1) ELISA Kit	1.56--100ng/mL	0.94ng/mL	Sandwich
E-EL-H2168	Human VWF (Von Willebrand Factor) ELISA Kit	1.56--100ng/mL	0.94ng/mL	Sandwich
E-EL-H0060	Human EGFR (Epidermal Growth Factor Receptor) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-H1593	Human TRAIL/TNFSF10 (Tumor Necrosis Factor Related Apoptosis Inducing Ligand) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-H2306	Human TNF- β (Tumor Necrosis Factor Beta) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-H0059	Human EGF (Epidermal Growth Factor) ELISA Kit	3.91--250pg/mL	2.34pg/mL	Sandwich
E-EL-H0764	Human F5 (Coagulation Factor V/Plasma Factor V) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-H0108	Human IFN- γ (Interferon Gamma) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-H0149	Human IL-1 β (Interleukin 1 Beta) ELISA Kit	7.81--500pg/mL	4.69pg/mL	Sandwich
E-EL-H0103	Human IL-10 (Interleukin 10) ELISA Kit	7.81--500pg/mL	4.69pg/mL	Sandwich
E-EL-H1277	Human HIF-1 α (Hypoxia Inducible Factor 1 Alpha) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-H1341	Human OPG (Osteoprotegerin) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-H0483	Human bFGF/FGF2 (Basic Fibroblast Growth Factor) ELISA Kit	12.5--800pg/mL	7.5pg/mL	Sandwich
E-EL-H0425	Human IGF1R (Insulin Like Growth Factor 1 Receptor) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-H0084	Human HGF (Hepatocyte Growth Factor) ELISA Kit	125--8000pg/mL	75pg/mL	Sandwich
E-EL-H0009	Human BAFF/CD257 (B-Cell Activating Factor) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-H0435	Human IGF-2 (Insulin Like Growth Factor 2) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-H0087	Human IGFBP-3 (Insulin Like Growth Factor Binding Protein 3) ELISA Kit	0.78--50ng/mL	0.47ng/mL	Sandwich
E-EL-H0048	Human IL-8 (Interleukin 8) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-H0068	Human FASL/TNFSF6 (Factor Related Apoptosis Ligand) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-H1629	Human EGFR2 (Epidermal Growth Factor Receptor 2) ELISA Kit	125--8000pg/mL	75pg/mL	Sandwich
E-EL-H2400	Human RANK (Receptor Activator of Nuclear Factor Kappa B) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-H0095	Human TNFSF14 (Tumor Necrosis Factor Ligand Superfamily, Member 14) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-H1603	Human VEGFR-2/KDR (Vascular Endothelial Growth Factor Receptor 2) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich

Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-H1133	Human ER α (Estrogen Receptor Alpha) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-H0067	Human FAS/CD95 (Factor Related Apoptosis) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-H1386	Human NF- κ B p105 (Nuclear Factor Kappa B p105) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-H0079	Human G-CSF (Granulocyte Colony Stimulating Factor 3) ELISA Kit	39.06--2500pg/mL	23.44pg/mL	Sandwich
E-EL-H0081	Human GM-CSF (Granulocyte Macrophage Colony Stimulating Factor) ELISA Kit	7.81--500pg/mL	4.69pg/mL	Sandwich
E-EL-H1846	Human PTGS2/COX-2 (Prostaglandin Endoperoxide Synthase 2) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-H1586	Human TGF- α (Transforming Growth Factor α) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-H2305	Human TACE (TNF α Converting Enzyme) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-H2132	Human STAT3 (Signal Transducer and Activator of Transcription 3) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-H2104	Human PAI1 (Plasminogen Activator Inhibitor 1) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-H2481	Human TNFSF13(Tumor Necrosis Factor Ligand Superfamily, Member 13)ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich
E-EL-H0010	Human BDNF (Brain Derived Neurotrophic Factor) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-H0666	Human CTNN β 1 (Catenin, Beta 1) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-H1451	Human MMP-9 (Matrix Metalloproteinase 9) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-M0049	Mouse TNF- α (Tumor Necrosis Factor Alpha) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-M0010	Mouse TNFRSF1A (Tumor Necrosis Factor Receptor Superfamily, Member 1A) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0643	Mouse TNFRSF1B (Tumor Necrosis Factor Receptor Superfamily, Member 1B) ELISA Kit	7.81--500pg/mL	4.69pg/mL	Sandwich
E-EL-M0044	Mouse IL-6 (Interleukin 6) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-M2441	Mouse EGFR (Epidermal Growth Factor Receptor) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-M0034	Mouse IGF-1 (Insulin-like growth factor 1) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich
E-EL-M0051	Mouse TGF- β 1 (Transforming Growth Factor β 1) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-M1215	Mouse TP53 (Tumor Protein 53) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-M0037	Mouse IL-1 β (Interleukin 1 Beta) ELISA Kit	7.81--500pg/mL	4.69pg/mL	Sandwich
E-EL-M0306	Mouse FV (Coagulation Factor V) ELISA Kit	3.13--200ng/mL	1.88ng/mL	Sandwich
E-EL-M1084	Mouse TRAIL/TNFSF10 (Tumor Necrosis Factor Related Apoptosis Inducing Ligand) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0048	Mouse IFN- γ (Interferon Gamma) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M1247	Mouse vWF (von Willebrand Factor) ELISA Kit	1.25-80ng/mL	0.75ng/mL	Sandwich
E-EL-M0046	Mouse IL-10 (Interleukin 10) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0025	Mouse EGF (Epidermal growth factor) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0170	Mouse bFGF/FGF2 (Basic Fibroblast Growth Factor) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich

Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-M0644	Mouse RANKL (Receptor Activator of Nuclear Factor Kappa B Ligand) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0687	Mouse HIF-1 α (Hypoxia Inducible Factor 1 Alpha) ELISA Kit	125--8000pg/mL	75pg/mL	Sandwich
E-EL-M0309	Mouse FX (coagulation factor X) ELISA Kit	0.31--20 μ g/mL	0.19 μ g/mL	Sandwich
E-EL-M0081	Mouse OPG (Osteoprotegerin) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-M0704	Mouse IGF-2 (Insulin Like Growth Factor 2) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-M2614	Mouse INS(Insulin) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-M0032	Mouse GM-CSF (Granulocyte-Macrophage Colony Stimulating Factor) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0031	Mouse G-CSF (Granulocyte Colony Stimulating Factor 3) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0484	Mouse FAS/CD95 (Factor Related Apoptosis) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-M0035	Mouse IGFBP-3 (Insulin-like Growth Factor Binding Protein 3) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-M0836	Mouse NF- κ B p105 (Nuclear factor NF-kappa-B p105 subunit) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-M1190	Mouse TGF- α (Transforming Growth Factor α) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0664	Mouse HGF (Hepatocyte Growth Factor) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-M0959	Mouse PTGS2/COX-2 (Prostaglandin Endoperoxide Synthase 2) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-M0028	Mouse FASL (Factor Related Apoptosis Ligand) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-M0476	Mouse ER α (Estrogen Receptor Alpha) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-M0195	Mouse BMP-6 (Bone Morphogenetic Protein 6) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich
E-EL-M0649	Mouse VEGFR-2/KDR (Vascuoar Endothelial Growth Factor Receptor 2) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-M2445	Mouse MCSF (Macrophage Colony Stimulating Factor 1) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0653	Mouse RANK (Receptor Activator of Nuclear Factor Kappa B) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-M0006	Mouse MCP-1 (Monocyte Chemotactic Protein 1) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0914	Mouse PAI-1 (Plasminogen Activator Inhibitor 1) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-M0117	Mouse CDKN2A (Cyclin Dependent Kinase Inhibitor 2A) ELISA Kit	39.06--2500pg/mL	23.44pg/mL	Sandwich
E-EL-M0036	Mouse IL-1 α (Interleukin 1 Alpha)ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0627	Mouse MMP-9 (Matrix Metalloproteinase-9) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich
E-EL-M0238	Mouse CASP3 (Caspase 3) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-M0042	Mouse IL-2 (Interleukin 2) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0053	Mouse CRP (C-Reactive Protein) ELISA Kit	0.39--25ng/mL	0.23ng/mL	Sandwich
E-EL-M0175	Mouse Bcl-2 (B-Cell Leukemia/Lymphoma 2) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0248	Mouse CTNNb1 (Catenin, Beta 1) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich

Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-M0298	Mouse CD40 (Cluster of Differentiation 40) ELISA Kit	0.78--50ng/mL	0.47ng/mL	Sandwich
E-EL-M0307	Mouse FVIII (Coagulation Factor VIII) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-M0815	Mouse NGF (Nerve growth factor) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-M0636	Mouse SCF (Stem Cell Factor) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0039	Mouse LEP (Leptin) ELISA Kit	125--8000pg/mL	75pg/mL	Sandwich
E-EL-M0203	Mouse BDNF (Brain Derived Neurotrophic Factor) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-M1057	Mouse STAT3 (Signal Transducer And Activator Of Transcription 3) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-R0019	Rat TNF- α (Tumor Necrosis Factor Alpha) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich
E-EL-R1091	Rat TNFRSF1A (Tumor Necrosis Factor Receptor Superfamily, Member 1A) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R1036	Rat TNFRSF1B (Tumor Necrosis Factor Receptor Superfamily, Member 1B) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0015	Rat IL-6 (Interleukin 6) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-R2386	Rat EGFR (Epidermal Growth Factor Receptor) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich
E-EL-R0010	Rat IGF-1 (Insulin-Like Growth Factor 1) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-R0084	Rat TGF- β 1 (Transforming Growth Factor β 1) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-R0012	Rat IL-1 β (Interleukin 1 Beta) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-R0230	Rat FV (Coagulation Factor V) ELISA Kit	3.13--200ng/mL	1.88ng/mL	Sandwich
E-EL-R1041	Rat TRAIL/TNFSF10 (Tumor Necrosis Factor Related Apoptosis Inducing Ligand) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0009	Rat IFN- γ (Interferon Gamma) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-R1079	Rat VWF (Von Willebrand Factor) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-R0016	Rat IL-10 (Interleukin 10) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-R0369	Rat EGF (Epidermal Growth Factor) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0091	Rat bFGF/FGF2 (Basic Fibroblast Growth Factor) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0841	Rat RANKL (Receptor Activator of Nuclear Factor Kappa B Ligand) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0513	Rat HIF-1 α (Hypoxia Inducible Factor 1 Alpha) ELISA Kit	7.81--500pg/mL	4.69pg/mL	Sandwich
E-EL-R0226	Rat F X (Coagulation Factor X) ELISA Kit	0.63--40ng/mL	0.38ng/mL	Sandwich
E-EL-R0050	Rat OPG (Osteoprotegerin) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-R0530	Rat IGF-2 (Insulin Like Growth Factor 2) ELISA Kit	46.88--3000pg/mL	28.16pg/mL	Sandwich
E-EL-R2466	Rat INS (Insulin) ELISA Kit	3.13--200ng/mL	1.88ng/mL	Sandwich
E-EL-R0008	Rat GM-CSF (Granulocyte-Macrophage Colony Stimulating Factor) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0454	Rat G-CSF (Granulocyte Colony Stimulating Factor) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0373	Rat FAS/CD95 (Factor Related Apoptosis) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich

Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-R0533	Rat IGFBP-3 (Insulin-Like Growth Factor Binding Protein 3) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich
E-EL-R0673	Rat NF- κ B p105 (Nuclear factor NF-kappa-B p105 subunit) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-R1012	Rat TGF- α (Transforming Growth Factor Alpha) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0496	Rat HGF (Hepatocyte Growth Factor) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-R0792	Rat PTGS2/COX-2 (Prostaglandin Endoperoxide Synthase 2) ELISA Kit	125--8000pg/mL	75pg/mL	Sandwich
E-EL-R2373	Rat FASL (Factor Related Apoptosis Ligand) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-R0118	Rat BMP-6 (Bone Morphogenetic Protein 6) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-R0865	Rat TACE (TNF α Converting Enzyme) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-R1060	Rat VEGFR-2/KDR (Vascuoar Endothelial Growth Factor Receptor 2) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-R0152	Rat RANK (Receptor Activator of Nuclear Factor Kappa B) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-R0633	Rat MCP-1 (Monocyte Chemotactic Protein 1) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-R0744	Rat PAI-1 (Plasminogen Activator Inhibitor 1) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-R0293	Rat CDKN2A (Cyclin Dependent Kinase Inhibitor 2A) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-R0011	Rat IL-1 α (Interleukin 1 Alpha) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0624	Rat MMP-9(Matrix Metalloproteinase 9)ELISA Kit	0.78--50ng/mL	0.47ng/mL	Sandwich
E-EL-R0160	Rat CASP3 (Caspase 3) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-R0013	Rat IL-2 (Interleukin 2) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-R0022	Rat CRP (C-Reactive Protein) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-R0096	Rat BCL-2 (B-Cell Leukemia/Lymphoma 2) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-R0170	Rat CTNNb (Catenin, Beta) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-R2377	Rat F VIII (Coagulation Factor VIII) ELISA Kit	7.81--500pg/mL	4.69pg/mL	Sandwich
E-EL-R0652	Rat NGF (Nerve Growth Factor) ELISA Kit	39.06--2500pg/mL	23.44pg/mL	Sandwich
E-EL-R0919	Rat SCF (Stem Cell Factor) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0582	Rat LEP (Leptin) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-R1235	Rat BDNF (Brain Derived Neurotrophic Factor) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-R1226	Rat STAT3 (Signal Transducer and Activator of Transcription 3) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich


Featured Publications

Cat. No.	Targets	Species	Publications
E-EL-H2532 E-EL-H0085	IFN- α (Interferon Alpha) IFN- β (Interferon Beta)	Human	Yang Q, Fu S, Wang J. Hepatitis C Virus Infection Decreases The Expression Of Toll-Like Receptors 3 And 7 Via Upregulation Of Mir-758[J]. Archives of Virology, 2014, 159(11): 2997-3003.
E-EL-H2443 E-EL-H2571 E-EL-H2338	IL-35(Interleukin 35) IL-37(Interleukin 37) IL-27(Interleukin 27)	Human	Gürkan A, Yücel AA, Sönmez C, et al. Serum Cytokine Profiles In Infants With Atopic Dermatitis[J]. Acta Dermatovenerol Croat, 2016, 24(4): 268-273.
E-EL-H2164 E-EL-H0483	VEGF-B(Vascular Endothelial Cell Growth Factor B) bFGF/FGF2(Basic Fibroblast Growth Factor)	Human	Ahmed F, Mohammed H. Effect Of Coronary Artery Bypass Grafting On The Plasma Levels Of VEGF, bFGF, Copper And NO As Angiogenic Biomarkers In Chronic Ischemic Heart Disease Patients[J]. Annals of Clinical and Laboratory Research, 2016, 4(3): 107-113.
E-EL-H1846	PTGS2/COX-2(Prostaglandin Endoperoxide Synthase 2)	Human	Aslan I, Özcan F, Karaarslan T, et al. Decreased Eicosapentaenoic Acid Levels In Acne Vulgaris Reveals The Presence Of A Proinflammatory State[J]. Prostaglandins & Other Lipid Mediators, 2016.
E-EL-H1779 E-EL-H1445 E-EL-H1451	PLAU/uPA(Urokinase-Type Plasminogen Activator) MMP-2(Matrix Metalloproteinase 2) MMP-9(Matrix Metalloproteinase 9)	Human	Chen H, Zou Y, Yang H, et al. Downregulation Of FoxM1 Inhibits Proliferation, Invasion And Angiogenesis Of Hela Cells In Vitro And In Vivo[J]. International Journal Of Oncology, 2014, 45: 2355-2364.
E-EL-H1615	sST2(soluble ST2)	Human	Yigitbasi M R, Guntas G, Atak T, et al. The Role Of Interleukin-33 As An Inflammatory Marker In Differential Diagnosis Of Idiopathic Granulomatous Mastitis And Breast Cancer[J]. Journal of Investigative Surgery, 2016.
E-EL-H1586	TGF- α (Transforming Growth Factor Alpha)	Human	Chen Y, Bai X, Qi Z, et al. The Hepatitis B Virus X Protein Promotes Pancreatic Cancer Through Modulation Of The PI3K/AKT Signaling Pathway[J]. Cancer Letters, 2016, 380(1): 98-105.
E-EL-H1451 E-EL-H2340	MMP-9(Matrix Metalloproteinase 9) 4-HNE(4-Hydroxynonenal)	Human	Wang D, Ma Y, Yang X, et al. Hypermethylation Of The Keap1 Gene Inactivates Its Function, Promotes Nrf2 Nuclear Accumulation, And Is Involved In Arsenite-Induced Human Keratinocyte Transformation[J]. Free Radical Biology and Medicine, 2015, 89: 209-219.
E-EL-H1451	MMP-9(Matrix Metalloproteinase 9)	Human	Mohamad S A, Gafar H S, Elhady M, et al. Cord Blood Matrix Metalloproteinase-9 As Predictor Of Respiratory Outcome In Preterm Infants[J]. Clinical Medicine and Diagnostics, 2015, 5(4): 53-59.
E-EL-H1451	MMP-9(Matrix Metalloproteinase 9)	Human	Yang X, Wang D, Ma Y, et al. Continuous Activation Of Nrf2 And Its Target Antioxidant Enzymes Leads To Arsenite-Induced Malignant Transformation Of Human Bronchial Epithelial Cells[J]. Toxicology & Applied Pharmacology, 2015, 289(2): 231-239.
E-EL-H1277	HIF-1 α (Hypoxia Inducible Factor 1 Alpha)	Human	Chandravanshi B, Bhonde R. Small Molecules Exert Anti-Apoptotic Effect And Reduce Oxidative Stressaugmenting Insulin Secretion In Stem Cells Engineered Isletsagainst Hypoxia[J]. European Journal of Pharmacology, 2016, 791: 424-432.
E-EL-H1277	HIF-1 α (Hypoxia Inducible Factor 1 Alpha)	Human	Chandravanshi B, Bhonde R R. Shielding Engineered Islets With Mesenchymal Stem Cells Enhance Survival Under Hypoxia[J]. Journal of Cellular Biochemistry, 2017.
E-EL-H0853 E-EL-H2578	CXCR1(CXC-Chemokine Receptor 1) CXCR2(CXC-Chemokine Receptor 2)	Human	Polat A, Tunc T, Erdem G, et al. Interleukin-8 And Its Receptors In Human Milk From Mothers Of Full-Term And Premature Infants[J]. Breastfeeding Medicine, 2016, 550(1): 346-356.

Cat. No.	Targets	Species	Publications
E-EL-H0649 E-EL-H0149 E-EL-H0753 E-EL-H0017	cTn-I/TNNI3(cardiac Troponin I) IL-1 β (Interleukin 1 Beta) NOS2/iNOS(Nitric Oxide Synthase 2, Inducible) CASP3(Caspase 3)	Human	Dwaich K H, Alamran F G, Alsheibani B I, et al. Melatonin Effects On Myocardial Ischemia-Reperfusion Injury: Impact On The Outcome In Patients Undergoing Coronary Artery Bypass Grafting Surgery[J]. International Journal of Cardiology, 2016, 221: 977-986.
E-EL-H0180	IL-9(Interleukin 9)	Human	Tuxun T, Apaer S, Ma H Z, et al. The Potential Role Of Th9 Cell Related Cytokine And Transcription Factors In Patients With Hepatic Alveolar Echinococcosis[J]. Journal of Immunology Research, 2015.
E-EL-H0119 E-EL-H0300 E-EL-H0002 E-EL-H0326 E-EL-H2164 E-EL-H1277	REN(Renin) AGT(Angiotensinogen) ACE(Angiotensin I Converting Enzyme) Ang-II (Angiotensin II) VEGF-B(Vascular Endothelial Cell Growth Factor B) HIF-1 α (Hypoxia Inducible Factor 1 Alpha)	Human	Nath M, Chandra P, Halder N, et al. Involvement Of Renin-Angiotensin System In Retinopathy Of Prematurity - A Possible Target For Therapeutic Intervention[J]. Plos One, 2016.
E-EL-H0111 E-EL-H1451	VEGF-A(Vascular Endothelial Cell Growth Factor A) MMP-9(Matrix Metalloproteinase 9)	Human	Zhang F, Lu Y, Chen Q, et al. Identification Of NCK1 As A Novel Downstream Effector Of STAT3 In Colorectal Cancer Metastasis And Angiogenesis[J]. Cellular Signalling, 2017.
E-EL-H0111	VEGF-A(Vascular Endothelial Cell Growth Factor A)	Human	He Y H, Li M F, Zhang X Y, et al. NLRC5 Promotes Cell Proliferation Via Regulating The AKT/VEGF-A Signaling Pathway In Hepatocellular Carcinoma[J]. Toxicology, 2016, 359: 47-57.
E-EL-H0109 E-EL-H0481 E-EL-H0661	TNF- α (Tumor Necrosis Factor Alpha) APAF1(Apoptosis Protease Activating Factor 1) CASP7(Caspase 7)	Human	Tang Y, Zhang Y C, Chen Y, et al. The Role Of miR-19b In The Inhibition Of Endothelial Cell Apoptosis And Its Relationship With Coronary Artery Disease[J]. Scientific Reports, 2015.
E-EL-H0109 E-EL-H0192	TNF- α (Tumor Necrosis Factor Alpha) IL-6R(Interleukin 6 Receptor)	Human	Veena V K, Popavath R N, Kennedy K, et al. In Vitro Antiproliferative, Pro-Apoptotic, Antimetastatic And Anti-Inflammatory Potential Of 2, 4-Diacetylphloroglucinol (DAPG) By Pseudomonas Aeruginosa Strain FP10[J]. Apoptosis, 2015, 20(10): 1281-1295.
E-EL-H0109 E-EL-H0103	TNF- α (Tumor Necrosis Factor Alpha) IL-10(Interleukin 10)	Human	Cheng P, Hua W, Zhang W, et al. Decreased Monocyte Subsets And TLR4-Mediated Functions In Patients With Acute Severe Fever With Thrombocytopenia Syndrome (SFTS)[J]. International Journal of Infectious Diseases, 2015, 43: 37-42.
E-EL-H0109 E-EL-H0102 E-EL-H0103	TNF- α (Tumor Necrosis Factor Alpha) IL-6(Interleukin 6) IL-10(Interleukin 10)	Human	Zhao Y, Jiang C, He J, et al. Multifunctional Dextran Sulfate-Coated Reconstituted High Density Lipoproteins Target Macrophages And Promote Beneficial Antiatherosclerotic Mechanisms[J]. Bioconjugate Chemistry, 2016.
E-EL-H0109 E-EL-H0102 E-EL-H0103	TNF- α (Tumor Necrosis Factor Alpha) IL-6(Interleukin 6) IL-10(Interleukin 10)	Human	Cheng P, Hua W, Zhang W J, et al. Inhibitory Effect Of miR-125b On Hepatitis C Virus Core Protein-Induced TLR2/MyD88 Signaling In THP-1 Macrophages[J]. World Journal of Gastroenterology, 2016, 22(17): 4354-4361.
E-EL-H0109 E-EL-H0150 E-EL-H0149	TNF- α (Tumor Necrosis Factor Alpha) IL-12(Interleukin 12) IL-1 β (Interleukin 1 Beta)	Human	Jiang Y, Kou J, Han X, et al. Ros-Dependent Activation Of Autophagy Through The PI3K/Akt/mTOR Pathway Is Induced By Hydroxysafflor Yellow A-Sonodynamic Therapy In THP-1 Macrophages[J]. Oxidative Medicine and Cellular Longevity, 2017.
E-EL-H0109 E-EL-H0102 E-EL-H0088	TNF- α (Tumor Necrosis Factor Alpha) IL-6(Interleukin 6) IL-1 α (Interleukin 1 Alpha)	Human	Chen W H, Luo G F, Lei Q, et al. Overcoming The Heat Endurance Of Tumor Cells By Interfering The Anaerobic Glycolysis Metabolism For Improved Photothermal Therapy[J]. Acs Nano, 2017, 11(2): 1419-1431.
E-EL-H0109	TNF- α (Tumor Necrosis Factor Alpha)	Human	Liu Q, Xu C, Ji G, et al. Sublethal Effects Of Zinc Oxide Nanoparticles On Male Reproductive Cells[J]. Toxicology in Vitro, 2016, 35: 131-138.
E-EL-H0109	TNF- α (Tumor Necrosis Factor Alpha)	Human	Zhang X P, Shang W L, Yuan J Z, et al. Positive Feedback Cycle Of TNF- α Promotes Staphylococcal Enterotoxin B-Induced THP-1 Cell Apoptosis[J]. Frontier in Cellular and Infection Microbiology, 2016.

Cat. No.	Targets	Species	Publications
E-EL-H0108 E-EL-H0109 E-EL-H0150 E-EL-H0102	IFN- γ (Interferon Gamma) TNF- α (Tumor Necrosis Factor Alpha) IL-12(Interleukin 12) IL-6(Interleukin 6)	Human	He W, Ren Y, Wang X, et al. C Reactive Protein And Enzymatically Modified LDL Cooperatively Promote Dendritic Cell-Mediated T Cell Activation[J]. Cardiovascular Pathology, 2017, 29: 1-6.
E-EL-H0108 E-EL-H0101 E-EL-H0105 E-EL-H0103	IFN- γ (Interferon Gamma) IL-4(Interleukin 4) IL-17(Interleukin 17) IL-10(Interleukin 10)	Human	Wang Z Y, Jiang M J, Wu M H, et al. Inflammatory Characteristics In Different Types Of Nonallergic Rhinitis[J]. International Journal of Clinical and Experimental Pathology, 2017, 10(3): 2887-2894.
E-EL-H0108	IFN- γ (Interferon Gamma)	Human	Wu R X, Yu Y, Yin Y, et al. Platelet Lysate Supports The In Vitro Expansion Of Human Periodontal Ligament Stem Cells For Cytotherapeutic Use[J]. Journal of Tissue Engineering & Regenerative Medicine, 2016.
E-EL-H0102 E-EL-H0109	IL-6(Interleukin 6) TNF- α (Tumor Necrosis Factor Alpha)	Human	Hao S Y, Ji J J, Zhao H T, et al. Mitochondrion-Targeted Peptide SS-31 Inhibited Oxidized Low-Density Lipoproteins-Induced Foam Cell Formation Through Both ROS Scavenging And Inhibition Of Cholesterol Influx In Raw264.7 Cells[J]. Molecules, 2015, 20(12): 21287-21297.
E-EL-H0102 E-EL-H0109	IL-6(Interleukin 6) TNF- α (Tumor Necrosis Factor Alpha)	Human	Liao Y, Liu P, Guo F, et al. Oxidative Burst Of Circulating Neutrophils Following Traumatic Brain Injury In Human[J]. Plos One, 2013.
E-EL-H0102 E-EL-H0048 E-EL-H0149	IL-6(Interleukin 6) IL-8(Interleukin 8) IL-1 β (Interleukin 1 Beta)	Human	Tian P, Li B, He C, et al. Antidiabetic (Type 2) Effects Of Lactobacillus G15 And Q14 In Rats Through Regulation Of Intestinal Permeability And Microbiota[J]. Food & Function, 2016.
E-EL-H0102	IL-6(Interleukin 6)	Human	Tang J N, Shen L, Liu C L, et al. Plasma Levels Of C1q/TNF-Related Protein 1 And Interleukin 6 In Patients With Acute Coronary Syndrome Or Stable Angina Pectoris[J]. American Journal of The Medical Sciences, 2015, 349(2): 130-136.
E-EL-H0101 E-EL-H0191 E-EL-H0104	IL-4(Interleukin 4) IL-5(Interleukin 5) IL-13(Interleukin 13)	Human	Li H T, Zhang T T, Chen Z G, et al. Intranasal Administration Of CpG Oligodeoxynucleotides Reduces Lower Airway Inflammation In A Murine Model Of Combined Allergic Rhinitis And Asthma Syndrome[J]. International Immunopharmacology, 2015, 28: 390-398.
E-EL-H0099 E-EL-H0101 E-EL-H0102 E-EL-H0103 E-EL-H0108	IL-2(Interleukin 2) IL-4(Interleukin 4) IL-6(Interleukin 6) IL-10(Interleukin 10) IFN- γ (Interferon Gamma)	Human	Li W, Dong H, Huang Y, et al. Clonorchis Sinensis Co-Infection Could Affect The Disease State And Treatment Response Of Hbv Patients[J]. Plos Neglected Tropical Diseases, 2016.
E-EL-H0099	IL-2(Interleukin 2)	Human	Zhou N, Fu Z, Li H, et al. Ketamine, As Adjuvant Analgesics For Patients With Refractory Cancer Pain, Does Affect IL-2/IFN- γ Expression Of T Cells In Vitro? A Prospective, Randomized, Double-Blind Study[J]. Medicine, 2017.
E-EL-H0099	IL-2(Interleukin 2)	Human	Zeng Z P, Wang K, Li Y Y, et al. Down-Regulation Of MicroRNA-451a Facilitates The Activation And Proliferation Of CD4+ T Cells By Targeting Myc In Patients With Dilated Cardiomyopathy[J]. The Journal of Biological Chemistry, 2016.
E-EL-H0089	IL-1ra/IL-1F3(Interleukin 1 Receptor Antagonist)	Human	Yao R, He Y, Zeng Z, et al. Protective Effect Of Adiponectin On Paraquat-Induced Pulmonary Fibrosis In Mice[J]. Molecular & Cellular Toxicology, 2015, 11(2): 247-255.
E-EL-H0088 E-EL-H0048	IL-1 α (Interleukin 1 Alpha) IL-8(Interleukin 8)	Human	Tang Y, Jin X, Xiang Y, et al. The lncRNA MALAT1 Protects The Endothelium Against ox-LDL-Induced Dysfunction Via Upregulating The Expression Of The miR-22-3p Target Genes CXCR2 And AKT[J]. Febs Letters, 2015, 589(20): 3189-3196.
E-EL-H0086	IGF-1(Insulin-like Growth Factor 1)	Human	Kini A, Fernandes C, Suryawanshi D. Effect Of Kifay On Performance, Insulin-Like Growth Factor-1, And Thyroid Hormones In Broilers[J]. Asian-Australasian Journal of Animal Sciences, 2016, 29(10): 1451-1457.
E-EL-H0085	IFN- β (Interferon Beta)	Human	Chen J, Xu W, Chen Y, et al. Matrix Metalloproteinase 9 Facilitates Hepatitis B Virus Replication Through Binding With IFNAR1 To Repress IFN/JAK/STAT Signaling[J]. Journal of Virology, 2017.

Cat. No.	Targets	Species	Publications
E-EL-H0085	IFN-β(Interferon Beta)	Human	Peng D, Wang Z, Huang A, et al. A Novel Function Of F-Box Protein FBXO17 In Negative Regulation Of Type I IFN Signaling By Recruiting PP2A For IFN Regulatory Factor 3 Deactivation[J]. Journal of Immunology, 2016.
E-EL-H0077 E-EL-H0102	Hepc(Hepcidin) IL-6(Interleukin 6)	Human	Li L, Zhang J, Zhao M, et al. Applied Value Of Monitoring Serum Hepcidin In Differential Diagnosis Of Infection Versus Tumor Fevers[J]. Journal of Huazhong University of Science and Technology-Medical Sciences, 2017, 37(2): 253-256.
E-AB-32609 E-EL-H0101 E-EL-H0103	PKC ζ Polyclonal Antibody IL-4(Interleukin 4) IL-10(Interleukin 10)	Human	Fan H, Li L, Zhang Y, et al. Pkcζ In Prostate Cancer Cells Represses The Recruitment And M2 Polarization Of Macrophages In The Prostate Cancer Microenvironment[J]. Tumor Biology, 2017.
E-EL-H0109	TNF-α(Tumor Necrosis Factor Alpha)	Human	Belvedere R, Bizzarro V, Parente L, et al. The Pharmaceutical Device Prisma® Skin Promotes In Vitro Angiogenesis Through Endothelial To Mesenchymal Transition During Skin Wound Healing[J]. International Journal of Molecular Sciences, 2017.
E-EL-H0483	bFGF/FGF2(Basic Fibroblast Growth Factor)	Human	Esnafoğlu E, Ayyıldız S N. Decreased Levels Of Serum Fibroblast Growth Factor-2 In Children With Autism Spectrum Disorder[J]. Psychiatry Research, 2017.
E-EL-H2557 E-EL-H5658 E-EL-H0016 E-EL-H0149	NLRP3(NACHT, LRR and PYD domains-containing protein 3) PYCARD(Apoptosis-associated speck-like protein containing a CARD) CASP1(Caspase 1) IL-1β(Interleukin 1 Beta)	Human	Isaza-Guzmán D M, Medina-Piedrahíta V M, Gutiérrez-Henao C, et al. Salivary Levels Of NLRP3 Inflammasome-Related Proteins As Potential Biomarkers Of Periodontal Clinical Status[J]. Journal of Periodontology, 2017.
E-EL-H0085	IFN-β(Interferon Beta)	Human	Qian W, Wei X, Guo K, et al. The C-Terminal Effector Domain Of Non-Structural Protein 1 Of Influenza A Virus Blocks IFN-β Production By Targeting TNF Receptor-Associated Factor 3[J]. Frontiers in Immunology, 2017.
E-EL-H2571	IL-37(Interleukin 37)	Human	Tawfik M G, Nasef S I, Omar H H, et al. Serum Interleukin-37: A New Player In Lupus Nephritis[J]. International Journal of Rheumatic Diseases, 2017.
E-EL-H1277	HIF-1α(Hypoxia Inducible Factor 1 Alpha)	Human	Mathew S A, Chandravanshi B, Bhonde R. Hypoxia Primed Placental Mesenchymal Stem Cells For Wound Healing[J]. Life Sciences, 2017.
E-EL-H0048 E-EL-H0109 E-EL-H2172 E-EL-H1640 E-EL-R0580 E-EL-R0504 E-EL-R0019 E-EL-R0488 E-EL-R0902	IL-8(Interleukin 8) TNF-α(Tumor Necrosis Factor Alpha) HO1(Heme Oxygenase 1) PECAM1/CD31(Platelet/Endothelial Cell Adhesion Molecule 1) LRP-1(Low-Density Lipoprotein-Receptor-Related Protein) HDL(High Density Lipoprotein) TNF-α(Tumor Necrosis Factor Alpha) HO1(Heme Oxygenase 1) sPECAM-1(Soluble Platelet Endothelial Cell Adhesion Molecule 1)	Human&Rat	Yan Z, Wang W, Wu Y, et al. Zinc Oxide Nanoparticle-Induced Atherosclerotic Alterations In Vitro And In Vivo[J]. International Journal of Nanomedicine, 2017, 12: 4433-4442.
E-EL-H0108 E-EL-H0101	IFN-γ(Interferon Gamma) IL-4(Interleukin 4)	Human	Xia G, Zheng X, Yao X, et al. Expression Of Programmed Cell Death-1 And Its Ligand B7 Homolog 1 In Peripheral Blood Lymphocytes From Patients With Peripartum Cardiomyopathy[J]. Clinical Cardiology, 2016.
E-EL-H0108 E-EL-H0101	IFN-γ(Interferon Gamma) IL-4(Interleukin 4)	Human	Xia G, Xin S, Zheng X, et al. Decreased Expression Of Programmed Death 1 On Peripheral Blood Lymphocytes Disrupts Immune Homeostasis In Peripartum Cardiomyopathy[J]. International Journal of Cardiology, 2016, 223: 842-847.



Universal Biotechnology Private Limited

Mob: +91 9910 73 8844 E-mail: info@unibiotech.in Web: www.unibiotech.in

Add: A-33, 2nd Floor, Ring Road, Rajouri Garden, New Delhi- 110027 (INDIA).