


Validated Antibodies for Canine and Feline Research




see more

GeneTex understands that veterinary researchers, as well as scientists who utilize canine and/or feline model systems, require high-quality antibodies to achieve accurate and reproducible results. For this reason, we are specifically testing an increasing number of our products for performance with dog and cat samples in addition to human, mouse, rat, and other animal species. All of the antibodies listed below have data supporting their use with dog and/or cat cells, cell lysates, or tissues.


For more information about these reagents, or if you are interested in testing these antibodies in previously untested species, please contact the GeneTex Technical Support Team at info@genetex.com.




Cat colon
ACE2 antibody [HL1092]
(GTX636265)



gamma Tubulin antibody
[HL1175] (GTX636480)




Dog brain
TLR3 antibody [HL2314]
(GTX638477)




MDCK cells
Vimentin antibody
(GTX100619)




Highly Cited Antibodies




ALDH1A1 antibody
(GTX123973)



Bcl-2 antibody [N1N2], N-term
(GTX100064)



Desmin antibody
(GTX103557)




ZO-1 antibody
(GTX108592)




Cat. No.	Product Name	Clonality	Reactivity	Applications
GTX11003	Actin antibody [AC-40]	Ms mAb	Bov, Carp, Chk, Dog, Goat, Gpig, Hm, Hu, Mk, Ms, Pig, Rat, Rb, Sheep, snail, Viper, Xenopus	ELISA, ICC/IF, IHC-Fr, WB
GTX11004	Actin antibody [AC-40] (Cy3)	Ms mAb	Bov, Carp, Chk, Dog, Goat, Gpig, Hm, Hu, Ms, Pig, Rat, Rb, Sheep, snail, Viper, Xenopus	ICC/IF, WB
GTX11005	Actin antibody [AC-40] (FITC)	Ms mAb	Bov, Carp, Chk, Dog, Goat, Gpig, Hm, Hu, Ms, Pig, Rat, Rb, Sheep, snail, Viper, Xenopus	ICC/IF, WB
GTX102419	Albumin antibody	Rb pAb	Cat, Dog, Hu, Ms, Rat	ICC/IF, IHC, IHC-Fr, IHC-P, WB
GTX123973	ALDH1A1 antibody	Rb pAb	Cat, Hu, Ms, Rat, Zfsh	ICC/IF, IHC-P, IHC-Wm, WB
GTX27291	alpha Tubulin antibody [DM1A]	Ms mAb	Bov, Chk, Dog, Fungi, Hu, Mk, Ms, Rat, Yeast	ICC/IF, IHC, IP, WB
GTX100064	Bcl-2 antibody [N1N2], N-term	Rb pAb	Cat, Dog, Hu, Ms, Rat, Zfsh	ICC/IF, IHC-P, WB
GTX109639	beta Actin antibody	Rb pAb	C. albicans, C. elegans, Cat, Chk, Dm, Dog, E. coli, Hm, Hrs, Hu, Mk, Mosq, Ms, Pig, Plant, Rat, Rice, Sheep, Treeshrew, Worm, Xenopus, Yeast	ICC/IF, IHC-Fr, IHC-P, IP, WB
GTX101435	beta Catenin antibody [N1N2-2], N-term	Rb pAb	Cat, Dog, Hu, Ms, Rat, Rb, Zfsh	ChIP assay, FACS, ICC/IF, IHC, IHC-Fr, IHC-P, IHC-Wm, IP, WB
GTX110723	Caspase 8 antibody	Rb pAb	Cat, Dog, Hu, Ms, Rat	ICC/IF, IHC-P, WB
GTX112888	Caspase 9 antibody [N2C3]	Rb pAb	Dog, Hu, Ms, Rat	IHC-P, IP, WB
GTX76185	CD9 antibody [MM2/57]	Ms mAb	Bov, Cat, Dog, Frt, Hrs, Hu, Lma, Mink, Ms, Pig, Rb, RMK	Blocking, FACS, IHC-Fr, IHC-P, WB
GTX76182	CD9 antibody [MM2/57] (Biotin)	Ms mAb	Bov, Cat, Dog, Frt, Hrs, Hu, Lma, Mink, Pig, Rb, RMK	FACS, WB
GTX76183	CD9 antibody [MM2/57] (FITC)	Ms mAb	Bov, Cat, Dog, Frt, Hrs, Hu, Lma, Mink, Pig, Rb, RMK	FACS, WB
GTX76186	CD9 antibody [MM2/57] (PE)	Ms mAb	Bov, Cat, Dog, Frt, Hrs, Hu, Lma, Mink, Pig, Rb, RMK	FACS, WB
GTX108624	Cyclin D1 antibody [N1C3]	Rb pAb	Cat, Dog, Hu, Rat	ELISA, ICC/IF, IHC-P, IP, Sandwich ELISA, WB
GTX108585	Cytochrome C antibody	Rb pAb	Dm, Dog, Hu, Ms, Rat	ICC/IF, IHC-Fr, IHC-P, WB
GTX103557	Desmin antibody	Rb pAb	Cat, Chk, Dog, Hu, Ms, Pig, Rat	FACS, ICC/IF, IHC-Fr, IHC-P, WB
GTX100443	E-Cadherin antibody	Rb pAb	Dog, Hu, Ms, Rat, Zfsh	ICC/IF, IHC-P, IHC-Wm, IP, PLA, WB
GTX134462	ERK1/2 antibody	Rb pAb	Cat, Dog, Hu, Ms, Rat	WB
GTX100118	GAPDH antibody	Rb pAb	Bacteria, Bov, C. elegans, Chk, Dm, Dog, E. coli, Fish, Hm, Hu, Insect, milkfish, Mk, Mosq, Ms, Nematode, Pig, Pika, Plant, Rat, Rb, Yeast, Zfsh	ICC/IF, IHC-P, IP, WB
GTX635013	GAPDH antibody [GT239] - VetSignal	Ms mAb	Cat, Dog	ICC/IF, IHC-P, WB
GTX20015	Heavy chain cardiac Myosin antibody [3-48]	Ms mAb	Bov, Dog, Hu, Ms, Pig, Rat, Rb	ELISA, ICC/IF, IHC-Fr, WB
GTX128839	Integrin beta 1 / CD29 antibody	Rb pAb	Dog, Hu, Ms, Rat	ICC/IF, IHC-Fr, IHC-P, IP, WB
GTX30773	NeuN antibody [2Q158]	Ms mAb	Crocodile, Dog, Hu, Ms, Rat, Rb	ICC/IF, IHC-Fr, IHC-P, WB
GTX23364	P Glycoprotein antibody [C219]	Ms mAb	Bov, Carp, Dog, Hu, Ms, Rat, Zfsh	ELISA, FACS, ICC/IF, IHC, IHC-Fr, IHC-P, IP, WB
GTX110720	p38 MAPK antibody [N1C3-2]	Rb pAb	Dog, Hu, Ms, Rat	ICC/IF, IHC-P, WB
GTX102425	p63 antibody [N2C1], Internal	Rb pAb	Cat, Dog, Hu, Ms, Rat	ICC/IF, IHC-Fr, IHC-P, IP, PLA, WB
GTX26401	pan Cytokeratin antibody [PCK-26]	Ms mAb	Bov, Carp, Cat, Chk, Dog, Goat, Gpig, Hm, Hu, Lzd, Ms, Pig, Rat, Rb, Sheep, Snk	Blocking, Dot, FACS, ICC/IF, IHC, IHC-Fr, IHC-P, WB
GTX75144	PCNA antibody [PC10] (FITC)	Ms mAb	Atlantic salmon, Cat, Chk, CIMK, Dog, Frt, Hm, Hrs, Hu, Ms, Rat, Rb, RMK, Sheep, Xenopus	FACS, IHC
GTX101025	PTEN antibody	Rb pAb	Cat, Dog, Hu, Ms, Pig, Rat	ICC/IF, IHC-P, WB
GTX75742	TLR4 antibody [HTA125]	Ms mAb	Cat, Dog, Gpig, Hu, Ms, Pig, RMK	Blocking, FACS, ICC/IF, IHC-Fr, IHC-P, IP, WB
GTX11210	TLR4 antibody [HTA125] (Biotin)	Ms mAb	Dog, Gpig, Hu, Pig, RMK	FACS, ICC/IF
GTX28378	TLR4 antibody [HTA125] (FITC)	Ms mAb	Dog, Gpig, Hu, Pig, RMK	FACS, ICC/IF
GTX75740	TLR4 antibody [HTA125] (Low endotoxin, azide free)	Ms mAb	Cat, Dog, Gpig, Hu, Ms, Pig, RMK	Blocking, FACS, Functional Assay, ICC/IF, IP, Neutralizing/Inhibition, WB
GTX11227	TLR4 antibody [HTA125] (PE)	Ms mAb	Dog, Gpig, Hu, Pig, RMK	FACS, ICC/IF
GTX70255	TSG101 antibody [4A10]	Rec Ms mAb	Dog, Hm, Hu, Mk, Ms, Pig, Rat, Zfsh	ELISA, EM, FACS, ICC/IF, IHC, IHC (Free Floating), IHC-P, IP, WB
GTX113030	VCP antibody [N1N2], N-term	Rb pAb	Dog, Hu, Ms, Rat, Zfsh	ICC/IF, IHC-P, IHC-Wm, WB
GTX100619	Vimentin antibody	Rb pAb	Cat, Dog, Hu, Ms, Rat	ELISA, ICC/IF, IHC, IHC-Fr, IHC-P, IP, Sandwich ELISA, WB
GTX108592	ZO-1 antibody	Rb pAb	Bov, Dog, Hu, Mk, Ms, Pig	ICC/IF, IHC-Fr, IHC-P, WB


Recombinant Antibodies




Integrin beta 1 / CD29 antibody [HL1255] (GTX636657)



p53 antibody [HL2242] (GTX638291)



PCNA antibody [HL1965] (GTX637858)



STAT3 antibody [HL1134] (GTX636400)




Cat. No.	Product Name	Clonality	Reactivity	Applications
GTX635479	53BP1 antibody [HL275]	Rec Rb mAb	Cat, Dog, Hu, Ms, Rat	ICC/IF, IHC-P, WB
GTX636265	ACE2 antibody [HL1092]	Rec Rb mAb	AGMK, Cat, Dog, Hu	ICC/IF, IHC-P, WB
GTX637047	ADAM17 antibody [HL1572]	Rec Rb mAb	Cat, Dog, Hu, Ms	IHC-P, WB
GTX636413	AKT1 antibody [HL1142]	Rec Rb mAb	Cat, Dog, Hu, Rat	ICC/IF, IHC-P, WB
GTX637940	Bcl-XL antibody [HL2039]	Rec Rb mAb	Cat, Dog, Hu	ICC/IF, IHC-P, WB
GTX637403	beta Catenin antibody [HL1753]	Rec Rb mAb	Cat, Dog, Hu, Ms	ICC/IF, IHC-P, WB
GTX637085	B-Raf antibody [HL1606]	Rec Rb mAb	Cat, Dog, Hu, Ms, Rat	ICC/IF, IHC-P, WB
GTX637077	Calnexin antibody [HL1598]	Rec Rb mAb	Cat, Dog, Hu, Ms, Rat	IHC-P, WB
GTX637068	CD28 antibody [HL1589]	Rec Rb mAb	Cat, Dog, Hu	ICC/IF, IHC-P, WB
GTX637231	CD44 antibody [HL1650]	Rec Rb mAb	Cat, Dog, Hu	ICC/IF, IHC-P, WB
GTX635979	clAP1 antibody [HL1045]	Rec Rb mAb	Cat, Dog, Hu, Ms, Rat	IHC-P, WB
GTX636468	CPT1A antibody [HL1164]	Rec Rb mAb	Cat, Dog, Hu, Ms	IHC-P, WB
GTX636580	Dicer antibody [HL1232]	Rec Rb mAb	Cat, Dog, Hu, Ms, Rat	ICC/IF, IHC-P, WB
GTX637063	Dopamine Receptor D2 antibody [HL1584]	Rec Rb mAb	Cat, Dog, Hu, Ms, Rat	ICC/IF, IHC-P, WB
GTX636577	E-Cadherin antibody [HL1229]	Rec Rb mAb	Cat, Dog, Hu, Ms, Rat	ICC/IF, IHC-P, WB
GTX637446	EPO antibody [HL1794]	Rec Rb mAb	Cat, Dog, Hu	ELISA, IHC-P, Sandwich ELISA, WB
GTX636884	FCER1G antibody [HL1418]	Rec Rb mAb	Cat, Dog, Hu, Ms, Rat	ICC/IF, IHC-P, WB
GTX637074	FGFR3 antibody [HL1595]	Rec Rb mAb	Cat, Dog, Hu	WB
GTX636487	GAD65 + GAD67 antibody [HL1181]	Rec Rb mAb	Cat, Dog, Hu, Ms, Rat, Zfsh	ICC/IF, IHC-P, WB
GTX636493	GAD65 antibody [HL1187]	Rec Rb mAb	Cat, Dog, Hu, Ms, Rat, Zfsh	ICC/IF, IHC-P, WB
GTX636280	GAD67 antibody [HL1096]	Rec Rb mAb	Cat, Dog, Hu, Ms, Rat, Zfsh	IHC-P, IHC-Wm, WB
GTX637921	Galectin 3 antibody [HL2020]	Rec Rb mAb	Cat, Dog, Hu, Ms	ICC/IF, IHC-P, WB
GTX636480	gamma Tubulin antibody [HL1175]	Rec Rb mAb	Cat, Dm, Dog, Hu, Ms, Rat	ICC/IF, WB
GTX637966	GAPDH antibody [HL2062]	Rec Rb mAb	Cat, Dog, Hu, Ms, Rat, Zfsh	ICC/IF, WB
GTX638336	Glutamine synthetase antibody [HL2283]	Rec Rb mAb	Cat, Dog, Hu, Ms, Rat	IHC-P, WB
GTX637432	Heme Oxygenase 1 antibody [HL1780]	Rec Rb mAb	Cat, Dog, Hu, Ms, Rat	IHC-P, WB
GTX637378	HGS antibody [HL1730]	Rec Rb mAb	Cat, Dog, Hu	ICC/IF, IHC-P, WB
GTX637963	HMGB1 antibody [HL2059]	Rec Rb mAb	Cat, Dog, Hu, Ms, Rat	ELISA, ICC/IF, IHC-P, WB
GTX637916	HPRT antibody [HL2015]	Rec Rb mAb	Cat, Dog, Hu, Ms, Rat, Zfsh	ICC/IF, IHC-P, WB
GTX637059	Hsp70 antibody [HL1580]	Rec Rb mAb	Cat, Dog, Hu, Ms, Rat	ICC/IF, IHC-P, WB
GTX636657	Integrin beta 1 / CD29 antibody [HL1255]	Rec Rb mAb	Cat, Dog, Hu, Ms	ICC/IF, IHC-P, WB
GTX637918	ISG15 antibody [HL2017]	Rec Rb mAb	Cat, Hu, Ms, Rat	ICC/IF, IHC-P, WB
GTX637379	KDM6A antibody [HL1731]	Rec Rb mAb	Cat, Dog, Hu	ICC/IF, IHC-P, WB
GTX637019	MCL1 antibody [HL1544]	Rec Rb mAb	Cat, Dog, Hu, Ms, Rat	IHC-P, WB
GTX636826	Mre11 antibody [HL1386]	Rec Rb mAb	Dog, Hu, Ms, Rat	ICC/IF, IHC-P, WB
GTX637955	MX1 antibody [HL2051]	Rec Rb mAb	Cat, Dog, Hu, Ms	IHC-P, WB
GTX637005	NPEPPS antibody [HL1531]	Rec Rb mAb	Cat, Dog, Hu, Ms, Rat, Zfsh	ICC/IF, IHC-P, WB
GTX635797	p38 MAPK antibody [HL1006]	Rec Rb mAb	Cat, Dog, Hu, Ms, Rat	WB
GTX638291	p53 antibody [HL2242]	Rec Rb mAb	Cat, Dog, Hu, Ms	ICC/IF, IHC-P, WB
GTX637858	PCNA antibody [HL1965]	Rec Rb mAb	Cat, Dog, Hu, Ms, Rat, Zfsh	IHC-P, WB
GTX636799	PDE6D antibody [HL1359]	Rec Rb mAb	Cat, Dog, Hu, Ms, Rat, Zfsh	ICC/IF, IHC-P, WB
GTX637103	PI3 kinase C2 alpha antibody [HL1612]	Rec Rb mAb	Cat, Dog, Hu, Ms, Rat	IHC-P, WB
GTX637912	PIP5K1A antibody [HL2011]	Rec Rb mAb	Cat, Dog, Hu, Rat	ICC/IF, IHC-P, WB
GTX637979	Plasminogen antibody [HL2071]	Rec Rb mAb	Cat, Dog, Hu, Ms, Rat, Zfsh	ICC/IF, IHC-P, WB
GTX637951	PODXL antibody [HL2047]	Rec Rb mAb	Cat, Hu, Ms, Rat	IHC-P, WB
GTX636919	PP2A alpha antibody [HL1447]	Rec Rb mAb	Cat, Dm, Dog, Hu, Ms, Rat	ICC/IF, IHC-P, WB



Recombinant Antibodies





AKT1 antibody [HL1142]
(GTX636413)



beta Catenin antibody [HL1753]
(GTX637403)



CD44 antibody [HL1650]
(GTX637231)



E-Cadherin antibody [HL1229]
(GTX636577)



Cat. No.	Product Name	Clonality	Reactivity	Applications
----------	--------------	-----------	------------	--------------

	GTX636971	RAB5A antibody [HL1497]	Rec Rb mAb	Cat, Dog, Hu, Ms, Rat	ICC/IF, IHC-P, IP, WB
	GTX638546	RNase H1 antibody [HL2343]	Rec Rb mAb	Dog, Hu, Ms	IHC-P, WB
	GTX638544	RNase L antibody [HL2341]	Rec Rb mAb	Dog, Hu, Ms, Rat	WB
	GTX635976	RPB2 antibody [HL1042]	Rec Rb mAb	Cat, Dog, Hu, Ms, Rat, Zfsh	IHC-P, IHC-Wm, WB
	GTX638494	RPL32 antibody [HL2331]	Rec Rb mAb	Dog, Hu, Ms, Rat, Zfsh	ICC/IF, IHC-P, WB
	GTX636579	RRM1 antibody [HL1231]	Rec Rb mAb	Cat, Dog, Hu, Ms, Rat	ICC/IF, IHC-P, WB
	GTX638300	SDHB antibody [HL2251]	Rec Rb mAb	Dog, Hu, Ms, Rat	IHC-P, WB
	GTX638555	SERCA2 ATPase antibody [HL2352]	Rec Rb mAb	Cat, Dog, Hu, Ms, Rat	IHC-P, WB
	GTX637889	SLC1A5 antibody [HL1988]	Rec Rb mAb	Cat, Hu, Ms	ICC/IF, IHC-P, WB
	GTX637255	SNAP91 antibody [HL1657]	Rec Rb mAb	Cat, Hu, Ms, Rat	IHC-P, WB
	GTX637250	SOD1 antibody [HL1652]	Rec Rb mAb	Cat, Dog, Hu, Ms, Rat	ICC/IF, IHC-P, WB
	GTX636297	Somatostatin antibody [HL1101]	Rec Rb mAb	Cat, Dog, Ms, Rat	ICC/IF, IHC-Fr, IHC-P, WB
	GTX636800	StAR antibody [HL1360]	Rec Rb mAb	Cat, Dog, Hu, Ms	ICC/IF, IHC-P, WB
	GTX636400	STAT3 antibody [HL1134]	Rec Rb mAb	Cat, Dog, Hu, Ms, Rat	IHC-P, WB
	GTX637282	TCEB1 antibody [HL1683]	Rec Rb mAb	Cat, Dm, Dog, Hu, Ms	IHC-P, WB
	GTX636635	Thymidylate synthase antibody [HL1236]	Rec Rb mAb	Cat, Dog, Hu, Ms, Rat	ICC/IF, IHC-P, WB
	GTX638477	TLR3 antibody [HL2314]	Rec Rb mAb	Cat, Dog, Hu, Ms, Rat	IHC-P, WB
	GTX637405	VEGFA antibody [HL1755]	Rec Rb mAb	Cat, Dog, Hu, Ms	ICC/IF, IHC-P, WB
	GTX636980	Vimentin antibody [HL1506]	Rec Rb mAb	Cat, Hu, Ms, Rat	ICC/IF, IHC-P, WB
	GTX637622	Vinculin antibody [HL1873]	Rec Rb mAb	Cat, Dog, Hu, Ms, Rat, Zfsh	IHC-P, WB
	GTX638372	VPS34 antibody [HL2305]	Rec Rb mAb	Dog, Hu, Ms, Rat	ICC/IF, IHC-P, WB
	GTX635821	VPS35 antibody [HL1017]	Rec Rb mAb	Cat, Dog, Hu, Ms, Rat, Zfsh	ICC/IF, IHC-P, IP, WB
	GTX637962	VWF antibody [HL2058]	Rec Rb mAb	Cat, Hu	IHC-P, WB
	GTX637900	Wnt2 antibody [HL1999]	Rec Rb mAb	Hu, Ms, Rat, Zfsh	WB
	GTX636914	ZC3H12A antibody [HL1442]	Rec Rb mAb	Cat, Hu	WB
	GTX636491	ZO-1 antibody [HL1185]	Rec Rb mAb	Cat, Dog, Hu, Ms, Rat	ICC/IF, WB

Recombinant

Citation Support

KO/KD Validation

Comparable Abs

Orthogonal Validation

Protein Overexpression



Scan For
More Info

www.unibiotech.in

Universal Biotechnology Private Limited

+91 9910 73 8844 info@unibiotech.in