

DNA LADDER


Best Fit for your Choice

GeneDireX[®]


UB UNIVERSAL
BIOTECHNOLOGY

Enquiry  +91 9910 73 8844


READY TO USE




1.5 % TAE agarose gel
100bp DNA Ladder
(100-3000bp)




1 % TAE agarose gel
1kb DNA Ladder
(250-10,000bp)



2% TAE agarose gel
50bp DNA Ladder
(50-1500bp)



0.7 % TAE agarose gel
XLarge DNA Ladder
(250-25,000bp)



1.5 % TBE agarose gel
OneMARK 100bp DNA Ladder
(100-3000bp)

www.unibiotech.in